

(19) World Intellectual Property Organization  
International Bureau



(43) International Publication Date  
5 June 2003 (05.06.2003)

PCT

(10) International Publication Number  
**WO 03/046126 A2**

(51) International Patent Classification<sup>7</sup>: C12N

(21) International Application Number: PCT/US02/34001

(22) International Filing Date: 23 October 2002 (23.10.2002)

(25) Filing Language: English

(26) Publication Language: English

(30) Priority Data:  
60/335,115 24 October 2001 (24.10.2001) US  
10/140,545 7 May 2002 (07.05.2002) US  
60/391,367 25 June 2002 (25.06.2002) US

(71) Applicant: AMERSHAM BIOSCIENCES (SV) CORP  
[US/US]; 928 East Arques Avenue, Sunnyvale, CA 94085  
(US).

(72) Inventors: SAMARTZIDOU, Hrisi; 938G La Mesa Ter-  
race, Sunnyvale, CA 94086 (US). HOUTS, Thomas; 8485  
Burchell Road, Gilroy, CA 95020 (US). YANG, Wen; 2025  
California Street, #30, Mountain View, CA 94040 (US).  
BUI, Son; 1243 Stellar Way, Milpitas, CA 95035 (US).  
HARKINS, Timothy; 884 Myrtle Street, San Jose, CA  
95126 (US).

(74) Agent: RONNING, Royal, N., Jr.; Amersham Bio-  
sciences Corp, 800 Centennial Avenue, Piscataway, NJ  
08855 (US).

(81) Designated States (*national*): AE, AG, AL, AM, AT, AU,  
AZ, BA, BB, BG, BR, BY, BZ, CA, CH, CN, CO, CR, CU,  
CZ, DE, DK, DM, DZ, EC, EE, ES, FI, GB, GD, GE, GH,  
GM, HR, HU, ID, IL, IN, IS, JP, KE, KG, KP, KR, KZ, LC,  
LK, LR, LS, LT, LU, LV, MA, MD, MG, MK, MN, MW,  
MX, MZ, NO, NZ, OM, PH, PL, PT, RO, RU, SD, SE, SG,  
SI, SK, SL, TJ, TM, TN, TR, TT, TZ, UA, UG, UZ, VN,  
YU, ZA, ZM, ZW.

(84) Designated States (*regional*): ARIPO patent (GH, GM,  
KE, LS, MW, MZ, SD, SL, SZ, TZ, UG, ZM, ZW),  
Eurasian patent (AM, AZ, BY, KG, KZ, MD, RU, TJ, TM),  
European patent (AT, BE, BG, CH, CY, CZ, DE, DK, EE,  
ES, FI, FR, GB, GR, IE, IT, LU, MC, NL, PT, SE, SK,  
TR), OAPI patent (BF, BJ, CF, CG, CI, CM, GA, GN, GQ,  
GW, ML, MR, NE, SN, TD, TG).

**Published:**

- without international search report and to be republished  
upon receipt of that report
- with sequence listing part of description published sepa-  
rately in electronic form and available upon request from  
the International Bureau

*For two-letter codes and other abbreviations, refer to the "Guid-  
ance Notes on Codes and Abbreviations" appearing at the begin-  
ning of each regular issue of the PCT Gazette.*

(54) Title: ARTIFICIAL GENES FOR USE AS CONTROLS IN GENE EXPRESSION ANALYSIS SYSTEMS

(57) Abstract: Method of producing universal controls for use in gene expression analysis systems such as macroarrays, real-time PCR, northern blots, SAGE and microarrays. The controls are generated either from near-random sequence of DNA, or from inter-genic or intronic regions of a genome. Twenty-three specific control sequences are also disclosed. Also presented are methods of using these controls, including as negative controls, positive controls, and as calibrators of a gene expression analysis system.



WO 03/046126 A2

ARTIFICIAL GENES FOR USE AS CONTROLS IN GENE EXPRESSION  
ANALYSIS SYSTEMS

CROSS-REFERENCE TO RELATED APPLICATIONS

5

This application is a continuation-in-part of United States patent application number 10/140,545, filed May 7, 2002, which claims priority to United States provisional patent application number 60/289,202, filed May 10 7, 2001, and 60/312,420, filed August 15, 2001. This application also claims priority to United States provisional patent application serial number 60/335,115, filed October 24, 2001, and 60/391,367, filed June 25, 2002, the disclosures of which are incorporated herein by 15 reference in their entireties.

REFERENCE TO SEQUENCE LISTING SUBMITTED ON COMPACT DISC

The present application includes a Sequence Listing filed on one CD-R disc, provided in duplicate, 20 containing a single file named pto\_PB0181.txt, having 56 kilobytes, last modified on October 21, 2002, and recorded on October 21, 2002. The Sequence Listing contained in said file on said disc is incorporated herein by reference in its 25 entirety.

BACKGROUND OF THE INVENTION

1. Field of the Invention:

30

The present invention relates to a method of using artificial genes as universal controls in gene expression analysis systems. More particularly, the present invention relates to a method of producing universal Controls for use

in gene expression analysis systems such as macroarrays, real-time PCR, northern blots, SAGE and microarrays, such as those provided in the Microarray ScoreCard system.

5    2.    Description of Related Art:

Gene expression profiling is an important biological approach used to better understand the molecular mechanisms that govern cellular function and growth.

10    Microarray analysis is one of the tools that can be applied to measure the relative expression levels of individual genes under different conditions. Microarray measurements often appear to be systematically biased, however, and the factors that contribute to this bias are many and ill-

15    defined (Bowtell, D.L., *Nature Genetics* 21, 25-32 (1999); Brown, P.P. and Botstein, D., *Nature Genetics* 21, 33-37 (1999)). Others have recommended the use of "spikes" of purified mRNA at known concentrations as controls in microarray experiments. Affymetrix includes several for use

20    with their GeneChip products. In the current state of the art, these selected genes are actual genes selected from very distantly related organisms. For example, the human chip (designed for use with human mRNA) includes control genes from bacterial and plant sources.

25            Each of the prior art controls consists of transcribed sequences of DNA from some source. As a result, that source cannot be the subject of a hybridization experiment using those controls due to the inherent hybridization of the controls to its source. In addition,

30    the lack of universal references consistent from experiment to experiment and from species to species greatly reduces the ability for scientists to compare data across labs, users, or time. What is needed, therefore, is a set of universal controls that do not hybridize with the DNA of any

source which may be the subject of an experiment. More desirably, there is a need for a universal control for gene expression analysis which do not hybridize with any known source.

5

#### SUMMARY OF THE INVENTION

Accordingly, this invention provides a process of producing universal controls that are useful in gene  
10 expression analysis systems designed for any species and which can be tested to insure lack of hybridization with mRNA from sources other than the control DNA itself.

The invention relates in a first embodiment to a process for producing at least one universal control for use  
15 in a gene expression analysis system. The process comprises selecting at least one non-transcribed (preferably intergenic, also intronic) region of genomic DNA from a known sequence, designing primer pairs for said at least one non-transcribed region and amplifying said at least one non-  
20 transcribed region of genomic DNA to generate corresponding double stranded DNA, then cloning said double stranded DNA using a vector to obtain additional double stranded DNA and formulating at least one control comprising said double stranded DNA.

25 The present invention relates in a second embodiment to a process of producing at least one universal control for use in a gene expression analysis system wherein testing of said at least one non-transcribed region to ensure lack of hybridization with mRNA from sources other  
30 than said at least one non-transcribed region of genomic DNA is performed.

The present invention in a third embodiment relates to said process further comprising purifying said DNA and mRNA, determining the concentrations thereof and

formulating at least one control comprising said DNA or of said mRNA at selected concentrations and ratios.

Another embodiment of the present invention is a universal control for use in a gene expression analysis system comprising a known amount of at least one DNA generated from at least one non-transcribed region of genomic DNA from a known sequence, or comprising a known amount of at least one mRNA generated from DNA generated from at least one non-transcribed region of genomic DNA from a known sequence. The present invention may optionally include generating mRNA complementary to said DNA and formulating at least one control comprising said mRNA, by optionally purifying said DNA and mRNA, determining the concentrations thereof and formulating at least one control comprising said DNA or of said mRNA at selected concentrations and ratios.

Another embodiment of the present invention is a universal control for use in a gene expression analysis system wherein a known amount of at least one DNA sequence generated from at least one non-transcribed region of genomic DNA from a known sequence, a known amount of at least one mRNA generated from DNA generated from at least one non-transcribed region of genomic DNA from a known sequence is included, and the aforementioned control wherein, said DNA and mRNA do not hybridize with any DNA or mRNA from a source other than the at least one non-transcribed region of genomic DNA.

The present invention, relates to a method of using said universal control, as a negative control in a gene expression analysis system by adding a known amount of said control containing a known amount of DNA, to a gene expression analysis system as a control sample and subjecting the sample to hybridization conditions in the

absence of complementary labeled mRNA and examining the control sample for the absence or presence of signal.

Further, said controls can be used in a gene expression analysis system by adding a known amount of a said control containing a known amount of DNA to a gene expression analysis system as a control sample and  
5       subjecting the sample to hybridization conditions, in the presence of a said control containing a known amount of labeled complementary mRNA, and measuring the signal values  
10       for the labeled mRNA and determining the expression level of the gene transcript based on the signal value of the labeled mRNA.

Additionally, said controls may be used as calibrators in a gene expression analysis system by adding a  
15       known amount of a said control containing known amounts of several DNA sequences to a gene expression analysis system as control samples and subjecting the samples to hybridization conditions in the presence of a said control containing known amounts of corresponding complementary  
20       labeled mRNAs, each mRNA being at a different concentration and measuring the signal values for the labeled mRNAs and constructing a dose-response or calibration curve based on the relationship between signal value and concentration of each mRNA.

Also, the present invention relates to a method of  
25       using said controls as calibrators for gene expression ratios in a two-color gene expression analysis system by adding a known amount of at least one of said controls containing a known amount of DNA to a two-color gene  
30       expression analysis system as control samples and subjecting the samples to hybridization conditions in the presence of a said control containing known amounts of two differently labeled corresponding complementary labeled mRNAs for each DNA sample present and measuring the ratio of the signal

values for the two differently labeled mRNAs and comparing the signal ratio to the ratio of concentrations of the two or more differently labeled mRNAs.

A further embodiment of the present invention is a process of producing controls that are useful in gene expression analysis systems designed for any species and which can be tested to insure lack of hybridization with mRNA from sources other than the synthetic sequences of DNA from which the control is produced.

One or more such controls can be produced by a process comprising synthesizing a near-random sequence of non-transcribed DNA, designing primer pairs for said at least one near random sequence and amplifying said non-transcribed DNA to generate corresponding double stranded DNA, then cloning said double stranded DNA using a vector to obtain additional double stranded DNA and formulating at least one control comprising said double stranded DNA.

The process can also be used to produce at least one control for use in a gene expression analysis system wherein testing of said sequence of non-transcribed synthetic DNA to ensure lack of hybridization with mRNA from sources other than said sequence of non-transcribed DNA is performed.

Additionally, mRNA complementary to said synthetic DNA can be generated and formulated to generate at least one control comprising said mRNA.

DNA and mRNA can be subsequently purified, the concentrations thereof determined, and one or more controls comprising said DNA or said mRNA at selected concentrations and ratios be formulated.

Another embodiment of the present invention is a control for use in a gene expression analysis system produced by the process comprises synthesizing a near-random sequence of DNA, designing primer pairs for said synthetic

DNA and amplifying said DNA to generate corresponding double stranded DNA, then cloning said double stranded DNA using a vector to obtain additional double stranded DNA and formulating at least one control comprising a known amount of at least one said double stranded DNA or a known amount of at least one mRNA generated from said DNA, and optionally, wherein, said DNA and mRNA do not hybridize with any DNA or mRNA from a source other than said DNA sequence of non-transcribed DNA.

10           The present invention, additionally, relates to a method of using said controls containing a known amount of DNA, as a negative control in a gene expression analysis system including adding a known amount of said control containing a known amount of DNA to a gene expression analysis system as a control sample, and subjecting the sample to hybridization conditions in the absence of complementary labeled mRNA and examining the control sample for the absence or presence of signal.

20           Further, said controls may be used in a gene expression analysis system wherein a known amount of a said control containing a known amount of DNA is added to a gene expression analysis system as a control sample and subjecting the sample to hybridization conditions in the presence of a said control containing a known amount of labeled complementary mRNA and measuring the signal values for the labeled mRNA and determining the expression level of the gene transcript based on the signal value of the labeled mRNA.

30           The present invention, also relates to a method of using said controls as calibrators in a gene expression analysis system including adding known amounts of a said control containing known amounts of several DNAs to a gene expression analysis system as control samples and subjecting the samples to hybridization conditions in the presence of a



said control containing known amounts of corresponding complementary labeled mRNAs, each mRNA being at a different concentration and measuring the signal values for the labeled mRNAs and constructing a dose-response or  
5 calibration curve based on the relationship between signal value and concentration of each mRNA.

The present invention, additionally, relates to a method of using said controls as calibrators for gene expression ratios in a two-color gene expression analysis  
10 system comprising adding a known amount of at least one of said controls containing a known amount of DNA to a two-color gene expression analysis system as control samples and subjecting the samples to hybridization conditions in the presence of a said control containing known amounts of two  
15 differently labeled corresponding complementary labeled mRNAs for each DNA sample present and measuring the ratio of the signal values for the two differently labeled mRNAs and comparing the signal ratio to the ratio of concentrations of the two or more differently labeled mRNAs.

20 Further embodiments and uses of the current invention will become apparent from a consideration of the ensuing description.

#### BRIEF DESCRIPTION OF THE DRAWINGS

25 The above and other objects and advantages of the present invention will be apparent upon consideration of the following detailed description taken in conjunction with the accompanying drawings, in which like characters refer to  
30 like parts throughout, and in which:

FIG. 1 shows representative results for the selection of universal controls that do not cross-hybridize with human RNA;

FIG. 2 shows representative results for the selection of universal controls that do not cross-hybridization with each other;

FIG. 3 represents a performance evaluation of the  
5 universal controls;

FIG. 4 shows a scatter plot of raw signals for the calibration and ratio controls from a two-color hybridization experiment;

FIG. 5 shows calibration curves based on the  
10 Calibration controls for a representative hybridization experiment;

FIG. 6 presents the control nucleotide sequence of DR1 (SEQ ID NO: 1);

FIG. 7 presents the control nucleotide sequence of  
15 DR2 (SEQ ID NO: 2);

FIG. 8 presents the control nucleotide sequence of DR3 (SEQ ID NO: 3);

FIG. 9 presents the control nucleotide sequence of DR4 (SEQ ID NO: 4);

FIG. 10 presents the control nucleotide sequence  
20 of DR5 (SEQ ID NO: 5);

FIG. 11 presents the control nucleotide sequence of DR6 (SEQ ID NO: 6);

FIG. 12 presents the control nucleotide sequence  
25 of DR7 (SEQ ID NO: 7);

FIG. 13 presents the control nucleotide sequence of DR8 (SEQ ID NO: 8);

FIG. 14 presents the control nucleotide sequence of DR9 (SEQ ID NO: 9);

FIG. 15 presents the control nucleotide sequence  
30 of DR10 (SEQ ID NO: 10);

FIG. 16 presents the control nucleotide sequence of RC1 (SEQ ID NO: 11);

FIG. 17 presents the control nucleotide sequence of RC2 (SEQ ID NO: 12);

FIG. 18 presents the control nucleotide sequence of RC3 (SEQ ID NO: 13);

5 FIG. 19 presents the control nucleotide sequence of RC4 (SEQ ID NO: 14);

FIG. 20 presents the control nucleotide sequence of RC5 (SEQ ID NO: 15);

10 FIG. 21 presents the control nucleotide sequence of RC6 (SEQ ID NO: 16);

FIG. 22 presents the control nucleotide sequence of RC7 (SEQ ID NO: 17);

FIG. 23 presents the control nucleotide sequence of RC8 (SEQ ID NO: 18);

15 FIG. 24 presents the control nucleotide sequence of Utility1 (SEQ ID NO: 19);

FIG. 25 presents the control nucleotide sequence of Utility2 (SEQ ID NO: 20);

20 FIG. 26 presents the control nucleotide sequence of Utility3 (SEQ ID NO: 21);

FIG. 27 presents the control nucleotide sequence of Negative1 (SEQ ID NO: 22);

FIG. 28 presents the control nucleotide sequence of Negative2 (SEQ ID NO: 23);

25 FIG. 29 presents the nucleotide sequence of DR1s used in a spike mix (SEQ ID NO: 24);

FIG. 30 presents the nucleotide sequence of DR2s used in a spike mix (SEQ ID NO: 25);

30 FIG. 31 presents the nucleotide sequence of DR3s used in a spike mix (SEQ ID NO: 26);

FIG. 32 presents the nucleotide sequence of DR4s used in a spike mix (SEQ ID NO: 27);

FIG. 33 presents the nucleotide sequence of DR5s used in a spike mix (SEQ ID NO: 28);

FIG. 34 presents the nucleotide sequence of DR6s used in a spike mix (SEQ ID NO: 29);

FIG. 35 presents the nucleotide sequence of DR7s used in a spike mix (SEQ ID NO: 30);

5 FIG. 36 presents the nucleotide sequence of DR8s used in a spike mix (SEQ ID NO: 31);

FIG. 37 presents the nucleotide sequence of DR9s used in a spike mix (SEQ ID NO: 32);

10 FIG. 38 presents the nucleotide sequence of DR10s used in a spike mix (SEQ ID NO: 33);

FIG. 39 presents the nucleotide sequence of RC1s used in a spike mix (SEQ ID NO: 34);

FIG. 40 presents the nucleotide sequence of RC2s used in a spike mix (SEQ ID NO: 35);

15 FIG. 41 presents the nucleotide sequence of RC3s used in a spike mix (SEQ ID NO: 36);

FIG. 42 presents the nucleotide sequence of RC4s used in a spike mix (SEQ ID NO: 37);

20 FIG. 43 presents the nucleotide sequence of RC5s used in a spike mix (SEQ ID NO: 38);

FIG. 44 presents the nucleotide sequence of RC6s used in a spike mix (SEQ ID NO: 39);

FIG. 45 presents the nucleotide sequence of RC7s used in a spike mix (SEQ ID NO: 40);

25 FIG. 46 presents the nucleotide sequence of RC8s used in a spike mix (SEQ ID NO: 41);

FIG. 47 presents the nucleotide sequence of Utility1s used in a spike mix (SEQ ID NO: 42);

30 FIG. 48 presents the nucleotide sequence of Utility2s used in a spike mix (SEQ ID NO: 43);

FIG. 49 presents the nucleotide sequence of Utility3s used in a spike mix (SEQ ID NO: 44);

FIG. 50 presents the nucleotide sequence of Negativels used in a spike mix (SEQ ID NO: 45); and

FIG. 51 presents the nucleotide sequence of Negative2s used in a spike mix (SEQ ID NO: 46).

#### DETAILED DESCRIPTION OF THE INVENTION

5

The present invention teaches universal Controls for use in gene expression analysis systems such as microarrays. Many have expressed interest in being able to obtain suitable genes and spikes as controls for inclusion in their arrays.

10 An advantage of the universal Controls of this invention is that a single set can be used with assay systems designed for any species, as these Controls will not be present unless intentionally added. This contrasts with the concept of using genes from "distantly related species." For example, an analysis system directed at detecting human gene expression might employ a *Bacillus subtilis* gene as control, which may not be present in a human genetic material. But this control might be present in bacterial genetic material (or at least, cross hybridize), thus it may not be a good control for an experiment on bacterial gene expression. The novel universal Controls presented here provide an advantage over the state of the art in that the same set of controls can be used without regard to the species for the test sample RNA.

20 The present invention employs the novel approaches of using either non-transcribed genomic sequences or totally random synthetic sequences as a template and generating both DNA and complementary "mRNA" from such sequences, for use as controls. The Controls could be devised *de novo* by designing near-random sequences and synthesizing them resulting in synthetic macromolecules as universal controls. Totally synthetic random DNA fragments are so designed that they do not cross-hybridize with each other or with RNA from

any biologically relevant species (meaning species whose DNA or RNA might be present in the gene expression analysis system). The cost of generating such large synthetic DNA molecules can be high. However, they only need to be  
5 generated a single time. Additionally, fragment size can be increased by ligating smaller synthetic fragments together by known methods. In this way, fragments large enough to be easily cloned can be created. Through cloning and PCR sufficient quantities of DNA for use as controls can be  
10 produced and mRNA can be generated by *in vitro* transcription for use in controls.

A simpler approach is to identify sequences from the intergenic or intronic regions (referred here as non-transcribed regions) of genomic DNA from an organism, and  
15 use these as a template for synthesis via PCR (polymerase chain reaction). Ideally, sequences of around 1000 bases (could range from 500 to 2000 bases) are selected based on computer searches of publicly accessible sequence data. The criteria for selection include:

- 20 1. The sequence must be from a non-transcribed region; and
2. The sequence must not have homology with or be predicted to hybridize with any known / published gene or expressed sequence tag (EST).

25 PCR primer pairs are designed for the selected sequence(s) and PCR is performed using genomic DNA (as a template) to generate PCR fragments (double strand DNA) corresponding to the non-transcribed sequence(s) as the control DNA. Additional control DNA can be cloned using a  
30 vector and standard techniques. Subsequently, standard techniques such as *in vitro* transcription are used to generate mRNA (complementary to the cDNA and containing a poly-A tail) as the control mRNA. Standard techniques are

used for purifying the Control DNA and Control mRNA products, and for estimating their concentrations.

Empirical testing is also performed to ensure lack of hybridization between the Control DNA on the array and  
5 other mRNAs, as well as with mRNA from important gene expression systems (e.g., human, mouse, *Arabidopsis*, etc.).

The above approaches were used to generate twenty-three universal control sequences from intergenic regions of the yeast *Saccharomyces cerevisiae* genome. Specifically,  
10 using yeast genome sequence data publicly available (<http://genome-www.stanford.edu/Saccharomyces/>), intergenic regions approximately 1 kb in size were identified. These sequences were BLAST'd and those showing no homology to other sequences were identified as candidates for artificial  
15 gene controls. Candidates were analyzed for GC-content and a subset with a GC-content of  $\geq 36\%$  was identified. Specific primer sequences have been identified and primers synthesized. PCR products amplified with the specific primers have been cloned directly into the pGEM<sup>TM</sup>-T Easy  
20 vector (Promega Corp., Madison, WI). Both array targets and templates for spike mRNA have been amplified from these clones using distinct and specific primers.

A greater number of intergenic regions have been cloned for testing. DNA samples from all the candidates  
25 were amplified, spotted on glass microarray slides and hybridized with mRNA samples from several species and each candidate spike mRNA, respectively, to identify those that do not cross-hybridize. First, they were screened for no cross-hybridization with RNA from different biological  
30 species. mRNA from human (eight tissues: skeletal muscle, spleen, liver, heart, kidney, brain, placenta and lung), mouse (six tissues: skeletal muscle, spleen, liver, heart, kidney and brain), rat (six tissues: skeletal muscle, spleen, liver, heart, kidney and brain), yeast (*S.*

cerevisiae) and bacteria (*E. coli* and two Archaea species), as well as total RNA from plant (*Arabidopsis*, Oil Palm) were tested against the control candidates. Candidates that did not cross-react with the RNA samples from the species tested  
5 were then selected for cross-hybridization with each other. The candidates were hybridized with each candidate mRNA independently.

From the candidate clones that exhibited specific hybridization, twenty-three were included into the final set  
10 of universal controls. FIG. 6 through FIG. 28 presents the nucleotide sequences of the twenty-three controls spotted on the microarray slides, while FIG. 29 through 51 presents the nucleotide sequences of the twenty-three controls that were transcribed and used in a spike mix, respectively. SEQ ID  
15 NO: 1 through SEQ ID NO: 23 present the nucleotide sequences of the twenty-three controls spotted on the microarray slides, while SEQ ID NO: 24 through SEQ ID NO: 46 present the nucleotide sequences of the controls that were transcribed and used in a spike mix.

20 These universal controls, when included in microarray experiments, perform as:

1. Negative controls: Control DNA included in the array, but for which no complementary artificial mRNA is spiked into the RNA sample, serves as a  
25 negative control;
2. Calibration controls: Several different Control DNA samples may be included in an array, and the complementary Control mRNA for each is included at a known concentration, each having a  
30 different concentration of mRNA. The signals from the array features corresponding to these Controls or Calibrators may be used to construct a "dose-response curve" or calibration curve to



estimate the relationship between signal and amount of mRNA from the sample;

3. Ratio controls: In two-color microarray gene expression studies, it is possible to include different, known, levels of Control mRNA complementary to Control DNA in the labeling reaction for each channel. The ratio of signals for the two dyes from a particular gene can be compared to the ratio of signals from the two dyes of the Control mRNA. This can serve as a test of the accuracy of the system for determining gene expression ratios.

4. Utility controls: These controls can be added into the sample preparation steps (such as RNA extraction and purification) for normalization of the biological samples and assessment of sample losses during preparation. Alternatively, they can be added to labeling reactions as additional calibrators or ratios.

Mixtures of several different Control mRNA species can be prepared (spike mixes) at known concentrations and ratios to simplify and standardize the experimental protocol while providing a comprehensive set of precision and accuracy information. Table 1 demonstrates one embodiment of this concept. The mRNA from the final set of clones have been pre-mixed at specific concentrations and ratios so they can serve as the various controls when hybridized to their corresponding control DNA spotted on the arrays. Ten calibrators (those included in the labeling reaction at a ratio of 1:1) spanning a dynamic range of 4.5 orders of magnitude are included as calibration controls. Eight ratio controls are included, at two expression levels (low and medium to high) and reversed with respect to the reference and test samples.

The universal controls as shown in Table 1 can be used as references for microarray validation and standardization across biological species and experimental platforms. These controls can be used to verify the accuracy and precision of gene expression ratios, and the sensitivity and dynamic range of the microarray system. Through the use of Calibration (standard) curves, these controls may allow reporting gene expression levels in consistent mass units, improving the comparisons of results across laboratories.

The following examples demonstrate how these Control DNA and Control mRNA were generated, and then used as universal controls in microarray gene expression experiments. They are representative of the many different types of experiments that could benefit from the use of these controls. The following examples are offered by way of illustration and not by way of limitation.

Table 1. Suggested Control mRNA spike mix composition for two-color gene expression ratio experiments.

Control Type	Control Name	Target Cy3: Cy5 Ratio	mRNA in the Spike Mix (pg/2 $\mu$ l of spike)	
			Cy3	Cy5
Calibration	DR1s	1:1	30 000	30 000
Calibration	DR2s	1:1	10 000	10 000
Calibration	DR3s	1:1	3 000	3 000
Calibration	DR4s	1:1	1 000	1 000
Calibration	DR5s	1:1	300	300
Calibration	DR6s	1:1	100	100
Calibration	DR7s	1:1	30	30
Calibration	DR8s	1:1	10	10
Calibration	DR9s	1:1	3	3
Calibration	DR10s	1:1	1	1
Ratio	RC1s	3:1 low	300	100
Ratio	RC2s	1:3 low	100	300
Ratio	RC3s	3:1 high	3 000	1 000
Ratio	RC4s	1:3 high	1 000	3 000

Ratio	RC5s	10:1 low	300	30
Ratio	RC6s	1:10 low	30	300
Ratio	RC7s	10:1 high	10 000	1 000
Ratio	RC8s	1:10 high	1 000	10 000
Utility	utility1s	User defined	User defined	User defined
Utility	Utility2s	User defined	User defined	User defined
Utility	Utility3s	User defined	User defined	User defined
Negative	Negative1s	NA	0	0
Negative	Negative2s	NA	0	0

Example 1. Generation of Artificial Controls  
from Intergenic Regions of *S. cerevisiae* Genome.

Using yeast genomic sequence data publicly available (<http://genome-www.stanford.edu/Saccharomyces/>), intergenic regions (YIRs) approximately 1 kb in size were identified. These sequences were BLAST'd and those showing no homology to other sequences were identified as candidates for artificial gene controls. Candidates were analyzed for GC-content and a subset with a GC-content of  $\geq 36\%$  was identified. Specific primer sequences have been identified and synthesized. PCR products amplified with the specific primers have been cloned directly into the pGEM<sup>TM</sup>-T Easy vector (Promega Corp., Madison, WI). Both array targets and templates for spike mRNA have been amplified from these clones using distinct and specific primers.

When used as DNA controls, the YIR sequences were amplified by PCR with specific primers, using 5 ng of cloned

template (plasmid DNA) and a primer concentration of 0.5 $\mu$ M in a 100  $\mu$ l reaction volume, and cycled as follows: 35 cycles of 94°C 20 sec., 52°C 20 sec., 72°C 2 min., followed by extension at 72°C for 5 min.

5 All YIR control mRNAs for the spike mix are generated by *in vitro* transcription. Templates for *in vitro* transcription (IVT) are generated by amplification with specific primers that are designed to introduce a T7 RNA polymerase promoter on the 5' end and a polyT (T21) tail on  
10 the 3' end of the PCR products. Run-off mRNA is produced using 1  $\mu$ l of these PCR products per reaction with the AmpliScribe system (Epicentre, Madison, WI). IVT products are purified using the RNeasy system (Qiagen Inc., Valencia, CA) and quantified by spectrophotometry.

15 Initially, fifty intergenic region sequences have been cloned for testing. DNA samples from all the candidates were amplified, spotted on glass microarray slides and hybridized with mRNA samples from several species and each candidate spike mRNA, respectively, to identify  
20 those that do not cross-hybridize. First, they were screened for no cross-hybridization with RNA from different biological species. mRNA from human (8 tissues: skeletal muscle, spleen, liver, heart, kidney, brain, placenta and lung), mouse (6 tissues: skeletal muscle, spleen, liver,  
25 heart, kidney and brain), rat (6 tissues: skeletal muscle, spleen, liver, heart, kidney and brain), yeast (*S. cerevisiae*) and bacteria (*E. coli* and two Archaea species), as well as total RNA from plant (*Arabidopsis*, Oil Palm) were tested against the control candidates.

30 Figure 1 shows the hybridization of candidates with human brain mRNA. The results indicated that two YIR clones, 33 and 62, hybridized with human brain RNA while the other candidates did not (since no appreciable signal is detected). Clones, such as 33 and 62, that exhibited such

template (plasmid DNA) and a primer concentration of 0.5 $\mu$ M in a 100  $\mu$ l reaction volume, and cycled as follows: 35 cycles of 94°C 20 sec., 52°C 20 sec., 72°C 2 min., followed by extension at 72°C for 5 min.

5           All YIR control mRNAs for the spike mix are generated by *in vitro* transcription. Templates for *in vitro* transcription (IVT) are generated by amplification with specific primers that are designed to introduce a T7 RNA polymerase promoter on the 5' end and a polyT (T21) tail on  
10 the 3' end of the PCR products. Run-off mRNA is produced using 1  $\mu$ l of these PCR products per reaction with the AmpliScribe system (Epicentre, Madison, WI). IVT products are purified using the RNeasy system (Qiagen Inc., Valencia, CA) and quantified by spectrophotometry.

15           Initially, fifty intergenic region sequences have been cloned for testing. DNA samples from all the candidates were amplified, spotted on glass microarray slides and hybridized with mRNA samples from several species and each candidate spike mRNA, respectively, to identify  
20 those that do not cross-hybridize. First, they were screened for no cross-hybridization with RNA from different biological species. mRNA from human (8 tissues: skeletal muscle, spleen, liver, heart, kidney, brain, placenta and lung), mouse (6 tissues: skeletal muscle, spleen, liver,  
25 heart, kidney and brain), rat (6 tissues: skeletal muscle, spleen, liver, heart, kidney and brain), yeast (*S. cerevisiae*) and bacteria (*E. coli* and two Archaea species), as well as total RNA from plant (*Arabidopsis*, Oil Palm) were tested against the control candidates.

30           Figure 1 shows the hybridization of candidates with human brain mRNA. The results indicated that two YIR clones, 33 and 62, hybridized with human brain RNA while the other candidates did not (since no appreciable signal is detected). Clones, such as 33 and 62, that exhibited such

cross-hybridization were removed from the set of candidates for universal controls.

Candidates that did not cross-react with the RNA samples from the species tested were then tested for cross-hybridization with each other. The candidates were  
5 hybridized with each candidate mRNA independently. In Figure 2 the labeled mRNA made from clone #50 was specifically hybridized against all other candidate clones. It hybridized only to its corresponding target DNA and can be included  
10 into the candidate set. However, clone #52 bound to the spot of clone #49 besides its own and therefore was not included in the candidate set.

From the candidate clones that exhibited specific hybridization, twenty-three are included into the final set  
15 of universal controls. FIG. 6 through FIG. 28 presents the nucleotide sequences of the twenty-three controls spotted on the microarray slides, while FIG. 29 through 51 presents the nucleotide sequences of the twenty-three controls as used in a spike mix, respectively. The sequences of these clones are  
20 further presented in the Sequence Listing, incorporated herein by reference in its entirety, as follows:

SEQ ID NOs: 1 - 23 (nt, control nucleotide sequences, including calibration controls 1 through 10, ratio  
25 controls 1 through 8, utility controls 1 through 3, and negative controls 1 and 2 respectively);  
SEQ ID NOs: 24 - 46 (nt, spike mix nucleotide sequences, including calibration controls 1 through 10, ratio controls 1 through 8, utility  
30 controls 1 through 3, and negative controls 1 and 2 respectively);

Upon confirmation of the exact structure, each of the above-described nucleic acids of confirmed structure is recognized to be immediately useful as a control.

5                    Example 2. Performance Evaluation of  
                      the Artificial Controls.

                  The universal controls (both the spike mixes and their corresponding spotting samples) have been evaluated  
10 for their performance in real microarray experiments and tested for the following.

                  Experimental design, including array design and the hybridization sample concentration were tested (Figure 3). Control samples were spotted in five replicates and  
15 hybridized with probes prepared with the spike mix only or the spike mixes with skeletal muscle mRNA. The same array image in Figures 3 is shown at two different gray scales for easy visualization of signals across the entire dynamic range.

20                    Universal utility, including hybridization of the spikes on pre-arrayed slides from various species were also tested. The controls showed no cross-hybridization on human, rat, mouse, *Arabidopsis*, Yeast and *E. coli* pre-arrayed slides from commercial sources (data not shown).

25                    Spike mix performance was tested, including ratio performance and Calibration curves (Figures 4 and 5). The mRNA from the final set of clones have been pre-mixed at specific concentrations and ratios (see Table 1 above) so they can serve as the various controls when hybridized to  
30 their corresponding control DNA spotted on the arrays. Ten calibrators (those included in the labeling reaction at a ratio of 1:1), spanning a dynamic range of 4.5 orders of magnitude, are included as calibration controls. Eight ratio controls are included, at two expression levels (low and

medium to high) and reversed with respect to the reference and test samples.

Figure 4 shows a scatter plot of raw signals for the calibration and ratio controls from a two-color hybridization experiment. The Calibrators are accurately and precisely clustered at the 45-degree line and the ratios at their expected target values at high (labeled 'H') and low (labeled 'L') levels of expression.

Figure 5 shows calibration curves based on the Calibration controls for a hybridization experiment. In this "standard curve", the Cy3 and Cy5 signals from the calibration controls are plotted as a function of the amount of mRNA in the spike mix. The error bars represent the 95% confidence intervals for the mean value. From such curves, attributes such as the limit of detection, the linear dynamic range and the signal saturation limit can be assessed. The application of the universal controls for the generation of standard curves can be the first step towards true quantitation of expression levels from microarray experiments.

The controls as shown in Table 1 can be used as references for microarray validation and standardization across biological species and experimental platforms. These controls can be used to verify the accuracy and precision of gene expression ratios, and the sensitivity and dynamic range of the microarray system. Through the use of Calibration (standard) curves, these controls may allow reporting gene expression levels in consistent mass units, improving the comparisons of results across laboratories

The above examples illustrate specific aspects of the present invention and are not intended to limit the scope thereof in any respect and should not be so construed.

Those skilled in the art having the benefit of the teachings of the present invention as set forth above, can



effect numerous modifications thereto. These modifications are to be construed as being encompassed within the scope of the present invention as set forth in the appended claims.

What is claimed is:

1. A control for use in a gene expression analysis system comprising:
  - 5 (a) a known amount of at least one DNA selected from the group consisting of
    - (i) SEQ ID Nos: 1 - 23;
    - (ii) a degenerate variant of the sequence set forth in (i); and
    - 10 (iii) a complement of the sequence set forth in (i) and (ii); or
  - (b) a known amount of at least one mRNA transcribed from the group consisting of
    - (i) SEQ ID NOS: 24 - 46;
    - 15 (ii) a degenerate variant of the sequence set forth in (i); and
    - (iii) a complement of the sequence set forth in (i) and (ii).
- 20 2. A method of using a control as a negative control in a gene expression analysis system comprising:
  - adding a known amount of said control DNA of claim 1, to a gene expression analysis system as a control sample;
  - 25 subjecting the sample to hybridization conditions in the absence of complementary labeled mRNA;
  - examining the control sample for the absence or presence of signal.
- 30 3. A method of using controls in a gene expression analysis system comprising:
  - adding a known amount of said control DNA of claim 1, to a gene expression analysis system as a control sample;

subjecting the sample to hybridization conditions  
in the presence of a known amount of labeled  
complementary mRNA of claim 1;

5       measuring the signal values for the labeled mRNA  
and determining the expression level of the DNA based  
on the measured signal value.

4. A method of using controls as calibrators in a gene  
expression analysis system comprising:

10       adding a known amount of a said control containing  
known amounts of several DNAs of claim 1, to a gene  
expression analysis system as control samples;

      subjecting the samples to hybridization conditions  
in the presence of a said control containing known  
15       amounts of corresponding complementary labeled mRNAs of  
claim 1, each mRNA being at a different concentration;

      measuring the signal values for the labeled mRNAs  
and constructing a dose-response or calibration curve  
based on the relationship between signal value and  
20       concentration of each mRNA.

5. A method of using controls as calibrators for gene  
expression ratios in a two-color gene expression  
analysis system comprising:

25       adding a known amount of at least one of said  
controls containing a known amount of DNA of claim 1,  
to a two-color gene expression analysis system as  
control samples;

      subjecting the samples to hybridization conditions  
30       in the presence of a said control containing known  
amounts of two differently labeled corresponding  
complementary labeled mRNAs of claim 1, for each DNA  
sample present;

measuring the ratio of the signal values for the two differently labeled mRNAs and comparing the signal ratio to the ratio of concentrations of the two or more differently labeled mRNAs.

1/51




FIG. 1

2/51




FIG. 2




FIG. 3

4/51




FIG. 4






FIG. 5

6/51

DR1

nt: SEQ ID NO: 1

tgttggtccaagaaaggagggttttgttcatcagaaaagaattcagaaaa	50
gcaaggaaacagtactatcgtttagaatgtagaatgatagggttgcttgct	100
aattctattatggcacgaatgatacacccatattttcaacaaaatcaata	150
cccactagcatcattgagccaactatttgtcaatgcaaccattaccggta	200
cttcatcctgatttaacgagtctacttttttatcacgtcaaaatttactt	250
gttttcctgtaaaccgaaataaaggcaaaaaagacctgggtgcaattac	300
gaataaatgtacaataatcatcctgtttgcatagtaaacttccagttaga	350
gtcacacaacgcaatgaattttgacagttttctgtgcatattctttggt	400
aaacgtaaagaacaggcaacttttgggtacaatggattctagcccatatgg	450
ttcatttctgggtgcattcgcaaagtcagtatttgtctagctgtgttttct	500
ggctgagagacattatgatgttattcattgttatggatatctctgtagct	550
catgctgcttattttctccctaaaaaagttttttctctcgaatacattctt	600
gaccatttcatagtgaattcttgtacttattttaaaccaaaaatggaag	650
tattcatacatccccctatcaaaaacactcaataagtttcgaattattcg	700
ttcgtctaaacagtgtccaataactcaaaggggtattcaagacggcacaaa	750
atcagcatcttcccttatccgtgttccagaaataccacgctaagggtttt	800
cctcctacaatccataaaatcattaaggaggcagcttgaaaaatcttgaa	850
attcaaaagagttatcttgggctaatacgaaattaacgataaccagagtag	900
aatattcaagatcacagctccaccttagtttcgagg	936

FIG. 6

7/51

DR2

nt: SEQ ID NO: 2

gcgaataacccaaaacgagactactttttaccattacaaccattttctttt	50
tccctattttctcactgggttgacagaaatcagtggtgctatcatcctaccat	100
atgcgctaaacttattgtctttctcctcctagagatgctgtattccatgc	150
atattctgaacgatgggttggtgtttttatcaagcaaggttaatcacatg	200
gcgtggcttgctccacacatcagtagaaaacgcataccgcagcggaatcc	250
ttaaataataagtgattttactgttcatcaactacaatcggactctttca	300
caattacccttcttgttttccacatttactgttaaataagggatgtaca	350
gaaggcttaggaaaacctgtgctgaatactggatggacactgcattccca	400
cagtgaacttttatagatacactgtcagttattttcgaactttcatcaa	450
gttgctgagtttttagtatccctttgccttagctatatgtttgaatgagca	500
aaatatttgcaatgtctctagctttcttgaaatattggtttatattgagg	550
gcttggtgaagatttcaaatttcactttgaaatactcaggagaaaaatcat	600
gctcttttgataatttggtgactaaacatacataaaacagtttaattttg	650
ggtggtaatggctgtgtgactagctatagaaagaaaaaattaaaaaaaa	700
aaaaaaaaatcaagtagttcctgcactgcgacgtccattatagcattatg	750
aattgggtccctgatttacgcatgcgataaactattttttagcgcagccgca	800
tattatccgagaataacttccgacataagaaaattcgcagaaaatagata	850
aaaaactgctcttggcattcttcacttcctctattacacactgtgtcata	900
ccacaatcatctcacagtatgtatttgatgtttatacatgctataacgt	950
aaaacaatgtagaatatatatctaaatacctcacgggttttagtttagtgc	1000

FIG. 7

8/51

DR3

nt: SEQ ID NO: 3

tgtgagaaattcactagcttcacctaaaagcaaaaaatctcaaaattccc	50
aatattcacatagtcctaaagtaccgatagcaaccaacatatataaacagt	100
agtatttttacgaagctgaattgcaagattagtgagaggagaataaaccgga	150
taattttttttggttacgttattgttaaaggctataatattaggtgaaac	200
agaatgtcctagaagtttttttctttcatgttaaatttattgattcttgc	250
gcttcagctttttataaaaacataagaactgtttcttcacgttaacttcttg	300
tgccacatatataatgatgtactagtaatatgggtactatttggcagatgat	350
atttgattttttattcaagacgggttactgtttctacgattgatattttcat	400
tcctggatatcatcttgccagatcacttacaatttaggccgcgcctgaat	450
tgaagagtacttcaatacgtagtgtactgtccaaactctcttccaaattt	500
ttaatattttagctgggggttgggtaacaagtgagcaagggaaaaagtgaac	550
attttaagaagaacaataaaaatagcaagagatggaatggtaatgcttggc	600
tctcgagaagagtagcataaaacgagacttggtttaaaacaggatatgaca	650
tacttcaattcagctttccctatcagccgctcgagcagttatataggtgt	700
gttgccggagtaatttggcggaggccaacagtggctaggcggcaacgcct	750
ggaacacgcgctttaaagttctggaagggttcgcgaattgagaactgctca	800
ggggcgaaatacaggggcggccttggcggcagggggaggcctctgtgaag	850
ttagttatataagacttgctgtcatcgtttttttgatcccggcaggaact	900
atcttttattctcatacatacgggtcaagaagtataattatacataacata	950
gggacacggttcaggcaattgtccatatccc	980

FIG. 8

9/51

DR4

nt: SEQ ID NO: 4

gatgtctgttttctttatgcaggatattaaatacaagtttgtgcttaaga	50
aatttccatgaagatatcaacattttattgtaggaattgcataaatagatg	100
aattatgtgccgctggacgtttatagatagcataagcacaatgacttaaa	150
ggttataatactcattgatatcactctgattataaaaatcgtaatatgcga	200
atagggtgaactaatcggaataaccacatcgacacttcaagcttcaattct	250
atttcaactgtagtgctgctagtgagaatacaaaaagtagcatacgtga	300
tgtgcaaaaaatgcgctacttatcacacaagtaccttgcgcaagaagggt	350
actctaaaccggggccatcgcatcaccagacggagatgtattctttatga	400
agcaataattggagggtgtatcaagttcgaaactgctgatgctatggattt	450
acatctttctttatgcacaaggcttgcttggtttctgagtagttagtttt	500
tagatttttgtcaagtctggggtaagtttaattcgagcaaaaattaacggca	550
cgttatttctaatacatatgttggttcataatattcttttacaagaggtttg	600
gaatgatgtcaccgatgttagaatgttaggagaatttcatgtgaatttta	650
gtccaagtgttgaagttctcttctgcagttagggcacgtacatggcaacg	700
atatcgtttttgatgtattaatcttagtagggcgttgagtttgatgttac	750
ttttctcagggtgatgaagcgtgatgacgatgacaaaaatgggttataata	800
gggcgcactatcatcatgcgtgattgatatttaaccaatgtcttgagtac	850
atcaactccagaaaatgggtcatttatatgcctagcatgtattatttgaga	900
cataaagttttatctcgagaccttgacgtataggaaa	937

FIG. 9

10/51

DR5

nt: SEQ ID NO: 5

ctgtagatttgcatggacgcacgtcgcccatagcctaaactttggcaatg	50
atactcgttatttcgtaatatcagtcctgcaaggtgctgtgatttctctat	100
tttatattgcctattatTTTTTcaaatgatttgagccgttttaaattgag	150
tatgcaatgagtcTTTTgaatcaaccgtaaggcagttccataaccactgc	200
cacgaatacgtttcactaccttgaagaatctctaattgtaggccgtattct	250
tcgcacttagttctgacgatgtagacatctcattatataagagcataagc	300
gcctgtttctagaatcatttcttcgtgaccagctttttgagttatttcg	350
cggatTTTTgaaacatttctcgagcttgacgtgaacatccttatatttca	400
tgacaaactcgatcattggaacatccctgcctcgatttttagagctagtat	450
caaatttcaatctctttgtgatggagccccgctcctatttcaaagagaa	500
gtttcttgatgcatatgttattgaagtctgattatagcaagtgcaatgt	550
cgtctcaattatttttaactatttttagccatacatgttagttatcctcaa	600
agagagcctccagactgggaagcagtggttgtcatttcaaataagtagat	650
ttcacagtttgatgattttcgaagccaggattcattgggctttgagtaa	700
agagaagccgcgtattacgaacagcttacgatattgtaaaatattccctt	750
attgtggtgccccaatggatacatgccagagaaatgtctgtgaaattgaa	800
caattacaatgacgagagcaagtaatccggcggccttgctctctttcac	850
tagtaccgtctatatctcttgagcgccaatatgcgaaagctttcacaagg	900
ttgatgttcatgggtattcggcgtcgatagcgaattgctta	940

FIG. 10

11/51

DR6

nt: SEQ ID NO: 6

ttagtttggaaacagcagtgtagataccgtccttggatagagcgctggaga	50
tagctggtctcaatctggtggagtaccatgggacaccagtgatgactcta	100
gtgacttgatcagcgggaataccagtcaacatagtggtgaaatcaccgta	150
gttgaaaacagcttcagcaatttcaactgggtaagtttcagttggatgag	200
cagcttggaaatatagtagttcagccaaatgagctctgatatctgagacg	250
tagacacctaattcgaccaggttaactctttcgtcagagggagataaagt	300
agtgggtggctggggcagcagcgacaccagcagcaatagcagcgacaccag	350
caacaattgaagttagtttgaccatcttttttcgattgaactttttagat	400
cttttttagtgaagatgtgagctcactcgaatgtaaataacaatgccaaat	450
tgtcggaaagagttaatcaaagctgctctatcttatatgccgttttttaat	500
aagcgacggacgaacagataaattggttgaatagctatttcactgctgata	550
tttctcttacttgggctcccctatcccatactcttcaccactacaaatat	600
gcagttgccctttcttcaacaatgctttttttatagatctcgtatacggg	650
tccgcgcctttgtactacctatatcttattatgatataacaggagcaca	700
ggaatgttcggtacagggatgataacctttaaggaagttttggcatgcct	750
tgacaacttcaattaatctttggccaagaaaatgaaccagaaatcaaatt	800
ttattctgtgccctctgaacgagggcaatatccaatgtttgacactaaac	850
ggttgtcaggagaaaaattgaatgtttcccaaatcagaaacattaaaatc	900
cctctatatgatcagaggagtcgtacctgttaggggtatgagcgaggaaac	950

FIG. 11

12/51

DR7

nt: SEQ ID NO: 7

aatgagttaccgtctgttacttttgggacggtttttgcactaagaacaga	50
cgagtttacggttatcctcaacaagcaagcaagtatttgctaacttagat	100
gccattccgaatcattactcatacgttactattgagagatgttttacaat	150
agatgagaagaatacaatgtccagagctcctgggtatgctagagtgcata	200
tccaggtcttattcgaatcatatcataccgtccatttcaacaatggtgaa	250
atgtgggtccacatatatcagaaatcttaacatttagtgaggagagccagt	300
agaaaaatgtgcgcaagcggaaagaagtcattcacagacacgtttaacaa	350
aacaccaccacagcagccttgtctcttgattctgatcagtttgccatcga	400
agaagcaaaattgtgggtgttatttttttcaaacaaaacttttttggcaac	450
agcagttttcttctggatatttgtactttatcatccaaccgatgaaagct	500
ggtttcctgtcaacctacatttaaattggcccgtacttcttcaaaaccgct	550
agataagcaaattaaccaacttttgagcgtcctaaattccccttggctc	600
agaagactcgttaatatgggaagttaaagtcctaccatataatcaaattg	650
gaagctttctgtgttcgaatggctatttctaaccgctgggctattaatcag	700
aggggaagtgaaatgaccgagacgtattatacgtcatgttgacatcaaca	750
atttaaggaaaaaaataaaaaaagcaatgaaaaagggtttttttaagtt	800
gaagacccttttcaaataatgttgccttgaattgtatctaccgtctcgt	850
ttcttctgctttaccggtttttttttgccttctttagatatgtcttttatg	900
cttgaaagggtccggctttaatgcattcatctaaacgtagtattcctat	950
ttgaactgctaccaatccaccatgactttact	982

FIG. 12



13/51

DR8

nt: SEQ ID NO: 8

gtcaaaactccacagccaagtccaagagtacgaaaaaaaaacttttcacaac	50
gagttcaaaaaatgtgggatgaagactcccgtttttcaccgcctagaaaa	100
cagcgtttgctgagaaaaaaaaataaatcatcgagaagaagtatgtcatca	150
taggatgttcccattgtaaggtgatgtgtaacatactcgaacaaagaatg	200
tatagagctgaatattttctccttttaaatttcaaagaaaatgagaaggaaa	250
atctcaaacagaaaacttcgttcttttttctcaagtaagcaaaagcttattg	300
agacaaagcggataactacgatattaataacgttgatgaagctcgaaca	350
aagtttagcgtcgggttatgcttgccatatataaagatatatttgccctacat	400
tttcgttgaacgtagaatgatttttgcttttaataaaattttttgttggttc	450
tttcagtgccttcttcaactttgatacgaaagcaagtgcattagtagacaaca	500
agaactggccacaactatactatactcattttttcttgcccggtgttttaa	550
atgttttcatccacagcatttgatgggatgattggaagtgagacgttcga	600
gaaaatccatatttttgagtcaagaattcagataatatactgagatgatta	650
ggatggctgggttctacaaaaacacaaatatccggctagcaatgatcac	700
tgagcaaattaaagcgttaactcactcattattgtagcttatgcgtttct	750
cctcctctctttttttcctcgaaccggagtggaagatccaataacgtaat	800
attactgatgttggttattaaagctggcaaaaataacatgaggcgtaaaac	850
cgactgcggtaagatgaggggtataaggtggagatcaggcgaacaagctg	900
ttcta	905

FIG. 13

14/51

DR9

nt: SEQ ID NO: 9

aacgatccaatgagcggtttcatgatgccattgtttaatcagagtgatgaa	50
aaagaaatatttgcgaccttttttcggttacattgatcgtgaaattttaat	100
caaagataatataaggacgtgagatatttatctttttacttgaaattaac	150
aatagaattgcgctaagcggaataagagctttcgtaaaccctttctatattg	200
caccattgcgtaacgtataaaatggatgacctttacacaaacgcatgc	250
ttataatcttatgtttttcatagggtgtaatttggttgatgacgtagtct	300
aaatttgatgctatctgcaattgaggtacatataagaggtcaatttcggg	350
accaacccttttaatcgaaaaaacgtaattcactagggcaaggggagaac	400
ttagcagctaataatcgtaaaccctttcataactaaaaaatgcacttaccat	450
caacaaaaaactcaggaccaattttccaagcttttctaggtgattgcctat	500
aacacaaaaagattcgctcatacatgagatttttacatgtaatagcaatt	550
tgttccgatcagttgaaggctcatcaacgcacggcaggtacatccacacct	600
atcacaaagcccttcaataattcacctacgtaaagttataccgaaacatg	650
caaaatccatgaaaaattctgtatgataacgatcatatccttttgattg	700
gtggtacgatgctcaaagatagttattgttgacctgaggcaaaagcgga	750
aatgaaaaatccagatggggccaaaagcagaagtattgtgtacaacaatt	800
gcttcagcagttttaccaaaccgtttcccagcaatcatcaaaagttgcttt	850
agccacatttccgcaagatatctttgtggctcaacgaagagggtattcc	900
aatgcaa	908

FIG. 14

15/51

DR10

nt: SEQ ID NO: 10

aactttcccccggtttaccacattgaagctgggtgtggaagatttatttga	50
agaaactaaaacgtaccctgtcatttcctgagtccttcaacttagtg	100
tgaaagccgaacaattataatcctcggtagacaacagatttattgtacta	150
aagttactcttcctgttatcttccttgattttactgttatagcaatgacc	200
caccgcaatcaggagagccgccgtatggaatagcataccaagtcatataaa	250
tcgtcaacctatttaacgggggttcagggttctttttcagcgtagtagccctt	300
taacaagcgctgacaaagttgacactcagagaaaattcaggatttattgt	350
aatccagctactcatccttagatccgcttgaggcatgggtttttttcacc	400
ttgagaggctattttgggtaagccaggaaggctgaaaaatcccaaaggga	450
cacagtaataagaaattgttggtgtgtatgatgcatttagaactcaaaa	500
gacgagtttctgaaaatgcttacaatactccataggtaacatgatttttt	550
tattaaaaaagtatactgttcctttgggtaaaaaattatgcaacccttgag	600
tgtccgatgaagataagactacgaaacaatttgcggtaaattttttctgc	650
tattgacatttacacatgctccaatccattaccctttccattctcgtaat	700
aaaacctcgaactgttattttcatattttacatctagacgggtatcggcctc	750
aacaactccaaacaaaagtaaatagaaaagagccagacctatcgcacccgg	800
gtagagccagaaaatattttaaactatagttgacgtattctacggctgtt	850
gtttaggacaatactttttccttcacaggcttcgaattacgcacatgcag	900
aactcctgt	909

FIG. 15

16/51

RC1

nt: SEQ ID NO: 11

gtgagacctccggaatTTTgacgctgcaagtcaatctacgggaaagaaga	50
aatTTTTtaaacctaatagcaaaaataagctTTTcttggaataaagatttt	100
cggcaataaaaggtaaatgcagccaaaaatcaaaatacttcagaagaagt	150
cgtagcgaggactgctaaggggaagcggatttgaagatcctttccagaac	200
aagaaggagccgaaagctgtcaggaactgttcctgattTTTTtaggaaaac	250
aattaataggtatctcgtctagagtagtatctcgagcttccagaagttgc	300
agataatcaaaatcattgttttatccctTTTTtagattacagcttagaa	350
gagtagagagcaagtttactgaaacgggttccttgtttacaataatattcc	400
taacaaactTTTtacgaattaggatgcagcatgattTTTTtatattgcttcac	450
ttcctaaagtatgaatTTTTtatccgtagtcgcaacaaaacagctactgg	500
aaatctgcagcttgTTtaaaaaccggtagtttccgaatactcctcgtcctt	550
gagttgtataccgTTtaacttcctaggggtgcatgtgtctggcccaattg	600
gccacaaaatctggTcctattgacggTTTTctTTtgattttcagcatct	650
tcctctaagaaggacagaaaattatgtaatatatgggagaaacggcctcc	700
caactgctaagtgtccccggcagcacgagtaagcaaaattcaggcaaact	750
attgcattaagaagccgtacataattcagcgtgatatgatgaaattttgt	800
taattgcaaattTTtagtacgatttggttgtagtggtgtttatgcaagt	850
aattattgaaccctaagtagttactgtctctTTtgctgtaattcgtgga	900
ttcacggccctccagcaacatggattgaa	929

FIG. 16

17/51

RC2

nt: SEQ ID NO: 12

agtagcatttatgacccaaaagcgtacttaaattagcagcaaaaaaat	50
taataacgaaactataaggaaaatacgagggtactgattatgagagtcccc	100
gtttctcatttttgagacatgatctgaacaaggctgaaaacagcaatctt	150
tttcgataacttttgcaaaaatttcaaacattgttggttgatgcagcca	200
atTTTTtatagggtacagagcttaatgctttacatgtgctttatTTTtcggt	250
actttccttaaagtgtctacattatctctcaggacttgaatgtcttcggc	300
tgaattactataaaatcttgagttttctctgaagtttaatcctaagacaa	350
tagtggtgagtgatgtagttcacgtgtgtgccactggtaataatagagat	400
aactatctcagttaagtttgaaaaggtaaaaaatagtttaagtagtcatt	450
ttttgcgacggtcattcttctctgatgcacgttcttttagactacctataa	500
acaccattcttacggaattataatggaaataaaacatcagtacgtgttgc	550
tgtcgggtgatagaggggtaacagaaccttaattgaaaaattagcacagt	600
cataattttattaacatgattgttttctgtgggaaataagaaatttcagca	650
ccagtaaaagacgagaaatatagggcacataaatgcgctcttactcgtat	700
gttccaggatgaaaatgttttagggcatcaagtattgccgaaagggaata	750
tgctttaacaccagaaaatccactgtatactcgttacgggtaaacaaagc	800
aaaacgcagtgcggtgataatgttttctaaaatctctgcacactgttgaaat	850
gcggctctgatacttttagcccttagtacctgacgggtgcctaaaatgagga	900

FIG. 17

18/51

RC3

nt: SEQ ID NO: 13

tagttggaggttggtgagtaccagattgcttacaaaagaatagcgagcca	50
acatttgctctgcctcaggcctcttggtgctgcttgaagactcatcttat	100
atggccttttgatgtcatgatttgcttcttgatattatgtgttgatatta	150
aacaaattgatttttttttttttgcgatagcaagcagataatgaaagaga	200
caaggacttggaacatccgataagactgcgccgatatcgatcttacagtc	250
cttccttggtgtcatgactttcggaagagcatcctcgtcgactggtagtt	300
tgctgtctgtcacgtgctgaagggtctgatacatTTTTTTaaagataaga	350
gacggggtttacccttcggaggactaagcgagatctccaagtaaagatct	400
cgcttatcaagaaagcagccaagtgtggaacgtccttttttttggtttca	450
aaaagatattcaacagttttacactgcagctttaattgcctcaaaaggata	500
tcatgaggtgatctagggtcagaagggaagattacagcatcttgagttg	550
aatcacatctgcaaaagggtggtattattgacgttgctcttccttaatgga	600
aactcatgggggtttggaaaggaggtgcggtaatctatttttttcgaacac	650
aaaacctaaccttgaaaagaaactgtccaatttcattgaacttacctcag	700
aacgggccggagtctttgctttcagtcataacatggtctaattttcttcgaa	750
aagcttcatttaattgttagactgtggttttacaaggaaaaaacagtgctc	800
tatactgaagcgatacccgagaactaattaccttggtgtgacgattcggctc	850
agcgaaacggacatggtaaaattgggaatttgaaagcaggcagcagcctt	900
gtacagcgacatgacgataggtttagaatcccatcacgtacgagttgaa	950
g	951

FIG. 18

19/51

RC4

nt: SEQ ID NO: 14

tcctagagtagcgattcccccttcgcgtattcttacatcttcgaagagaac	50
ttctggtgtaagtataataaatattatagctctatcgaatggtgcaatta	100
tttaccaaattctcaataggaatccataatactacatacgataactaatat	150
tctagtatTTTTtatacttattatttctTTTTtattacaccagcaatcggt	200
gcaaattatcttctgatagaatttctgagggtatcctaaacttatgccat	250
tttcttggactgtaaatacacttggatgttggtgcattagtcaataatcg	300
gttcttgttccaacgattacatgtaaatagaaggagaaataattatggta	350
aatcatgcggcgggtccttttggatgcagtatccatagtcactacataa	400
caatcttagtcaccttgtattgattcaccacataatcctgcagagccgc	450
tatgtccttaatactgcgcgataactctcctacccctgaattttgagagcg	500
ccatagcaaaccgataaagctggcacaattaaaggatatcgggtgttgtag	550
aattaggtgcctcctgcttttttttttctcctgctcttatatccgttata	600
tccgaatgatttttatcgcttggttaaaaaatactttcccgatatatata	650
tatagtctcccttttaaatttggttccggtaagtttttaacaccaataaat	700
gaaaagaaatgactacggatgaatatgagccgcgcattgaatcaggtt	750
atgtaagtatcagaacccctaattatgatgtcactcttacccttcgatgg	800
ctaagcggcgactgggatgccgggaaaagctctacaaatctactaaaaaa	850
gtcaaatatacagctgtaaacttcttctcctcgtctacatcatggtaacga	900
ttgttcaatctttacttcgtgtcttttttttttctatgtactttctatt	950
ccaacctatgtgaagactaaaattcaccttagtaaacgtaaagacaatga	1000
cgataggtgc	1010

FIG. 19

20/51

RC5

nt: SEQ ID NO: 15

acaattgcatctgtgcttggctctcaagaacgggtgtttggtgcatcaaaa	50
gttttcgactgcttattttggtcggaaatataaaaaactcgatcctcttatac	100
taagcagtatacattcttctttttgaaatgaatgtactccgtaatatctt	150
cttattttggcattttcatccttaacttttgcattggctctgaactagtcag	200
atagttgcccttttcagcaaacctcttattattgaaagcatggtgtacat	250
ccgttatactattatattataagaaattgggatgccaatTTTTTTgcttt	300
tgTTTTgctgttttcttcttttgcgaaaagtaattgcagatttaatag	350
caggatattataccgttggtaaaaacttaaggattttatgaacaatagctt	400
caagtacagcattcatagaaccaactactaaggatgaaactagtatgttt	450
ttgtcaaaaatattttcttgaccttgctgtaacatcaagatctgtttctct	500
aagatattaaagttgagtaaaaacaaagctgatatgagaaaaatacgtaa	550
ttgctccacataatacgtgggtcagacataaaggtagaataacttgataca	600
gaagagattattcgggtactcttgatggcgtgcttgaactggtgcctctta	650
acaaccggtaatatagtcagatgagtcactacgagtgtgtgtagtagcaa	700
gtgttttacctacgtggcagtaagagtagctctatggttgtgtaatagtg	750
gtgcttattcctaatactctgaagtctgaagcgggtacagttggtctggtc	800
tatatcatgggtcaaaggagcaaacatatcttctgaagtgaccgcaaatag	850
tactatgatgtggttggaatataacttaaaaggaaataaccacaaggaa	900
ttgcacccatgtacacagtttttcccggaattgggaaaccagta	945

FIG. 20



21/51

RC6

nt: SEQ ID NO: 16

gcataatgcccgctataaaccttattttttatatgggggtctggcgcttc	50
gggaaaagagaggaaaacttgtaactcaatatatctcgatacaacattac	100
gttttgtaaatttatcacaaaagccaaatgatgatattctctcttgcaagt	150
tatcgaacattgattggtaatttggttgaaaattgttaatttattgaata	200
tttcttttgcaaaagaaatagtcctcagcgaaagctgggtacaaaatttac	250
atcatgagtttacgggatttgtaaatatcgctttttgcataaaaatacttt	300
gccgtttcccacccttgcatattcacttactcccccttcatatactcta	350
tgtaatgatgattaagctttggccgctaagtcctctcaattagtgttgatt	400
ttggttttattcatatgattcttcttttagtgaagtattgatcaattacgt	450
gagtcagctttttgaaaaccccatgttggaaggaattaggaaattattttg	500
cttactacgaccactaatttaccgccatttctgggcctttttattgacta	550
ttttgaccatgtgctcgactagaagaacggcatcataatctgctggtaga	600
gttagtctataatgattggttgaaaataaaggcataagagatattccacct	650
aaaattcaagttattgactttattatcaggatcttagtatccttttttgg	700
taagtcataattcaatgaactaggtctcgcaaactttttgttcgaaaagcg	750
gtagtgcatagttatgctaactctggatatatggcataaaccgtacaaca	800
ctagcccatttttttggaagtagtgagggcagctagactgtatgatgaat	850
attcgctgcatactgagttttttggtccttttttttatgtggctggcct	900
tacgatatg	909

FIG. 21

22/51

RC7

nt: SEQ ID NO: 17

gtttagtgcacccagactcgaatcttaaataccactttacacaccta	50
aatTTTgtcctcacaaatgaagacaggattcaaaaccgattaatagtag	100
cagaaactaaaaaagtacgaatattagtaaaattcatgttcttgaatcga	150
gctactatctttgtcgggagggtaaacgattataactcaaaatgactgga	200
actggtgattattaatTTTtacgtttcctgtgccaataagcgggaagataa	250
gaggatagaagaaaagaaaggcggcacttggcgaactacaatggcgatta	300
tattcatggcgattatattcatacaaaggtaatggaggcctcgataatg	350
gacaatattgagaaaatccttatgcttacttctcttaataaaaaatagac	400
acagccatttattatgcgtaaaaaagattaccacttgtcttcgatgcgt	450
gctgctgccaatcaacctTTTgagcgggaacttcgagctcgcaatgcgtct	500
ggaatgttgctagagacagtcttggttatctgtgacatgtgtttcgttca	550
ggcgtgtgagcatcttcttggttcgatttcaaaattaccgccttgactcgt	600
gaaactggataattcggttggcgttttcatataagtcgtctgatggcgaaa	650
acttttcccttacttagcatacagcaaatatccccatttgacggattttt	700
gaaaaatgagcccgcctaaccagaatgaactgcattaccaagcatttatg	750
taaacgttccgccaccatctttggtaagggtatactattatgttctggatt	800
taaggttgattcacaatTTTcatcaccaaaatctggtggcatgcctagt	850
tgtctggtttcaggcaatttagccatcatagaaaagcatcctctgtcttg	900
agttgagaaaatgttactcatagagccaaacaaataaacctgg	944

FIG. 22

23/51

RC8

nt: SEQ ID NO: 18

ttcttgccgttctcatttcctgcacagtttctttgattatgtttgcagaa	50
gaatttcctttatcgtttagtctaaacaaagaattcgttgtaaagaatttg	100
agagcggatcttgcatTTTTTatTTatcatgcttatgtTTTTtctttgat	150
gtaagaagaagcaagtaagatatgtgaatatcttatcactaattcaaata	200
actaagagagctcacaacgacaatttggtgacagcatgcgaagcaaagagc	250
agtgataccagtatctttcatccagtaataacatacgaactgatgttatag	300
ttaaatgttacattttgagagacttcaacctctcgaaaccaagaggttgg	350
ttttaactctgggtgacttcaagaagggtgggtaccttttaciaaagcttga	400
gacgaagcaatagtcagtcctctgtataacaaggagaccacctcattttcc	450
agtaactcttgaggcatgtcggatgggtttgccttgaataaaccgcagtca	500
ttataatgaatggcctgtactttcaaaacagtctggaaacagaaatccat	550
tgctgagggtaccttttagtagcactttcgtagtgaaggtttaaggtag	600
ttcttatttactgcacaagagtttacatttaaccactctaatagtaactg	650
ttagagtggtttaactgttaggtgatctgttcattccatttttcgtgttg	700
tatctcaagatgagatagccttagcggttgctacatacataaatctaaacat	750
ataaacacctgtgtaactcgttaacgtctgggcttccatgcttctacat	800
ttagaatgatgtagaccattttattccaagaggataagcaccctctgtgat	850
tcaaaatgataataagtgttgacgacaagttactctcgcagaattgttgt	900
caa	903

FIG. 23

24/51

## Utility1

nt: SEQ ID NO: 19

tgcaatttaaagagcggtacctgtaaataagaaggaagaacggttatgttat	50
taatggacttttttagtgtcatcgaattttatgtaatatataagaaggtag	100
aataatttggcaggataatgtgttagcaaaggaggaaatcgaataccttt	150
aaaagagaaaaaatTTTTtagctgcttaaatttctgtgttataccacccg	200
atagattttgagttatgctttctaattgatctgactgccaacgttttctt	250
tatgccatctgaattgtcaggaacaaagaagaaaaagaaagtTTTTtaa	300
aaatctgtggtcgtgtgtgatgtacctttcctttacatgcattaatgcgc	350
tctgaaatgtggtacgatatccttacagagaatatattttctgtatatcg	400
tgcaatggtgaataacctatgaaggaaagtacccatcgctcaaggtaagc	450
attccaggaggggtcgccagaaacttaaactagtttttagcgacagatccga	500
aaattgatagagacattgaaaaaatcactactccgtccttttttagtgctt	550
tctcaatgcataattttggtgcacgactaaaaaattctagaacactatag	600
ttgcatttttttgggccggaagaagaaaaacgcgatgtaactttaatgtcaa	650
ataaagttttcacctagtaagcgcgatacaaaaaaacacagaaatagcc	700
ataggaaagtgaattttgtcagccgactaaaattaaggtttagcttacaaa	750
gcagcaaaaaatttgacatcgcacggtattccctgaaaaaggagcaggca	800
ggtgctgtatatTTTTTTTcggttcctgcctcttacatggcgtcgggtgtat	850
cttaaataactaaagtgagctgactacccttttgagtgccttatgtgacct	900
ctgatctcgaaagtaaacagagatacctaatttcacagccacttttggt	950
gcggacactgacgggatgtgttg	973

FIG. 24

25/51

Utility2

nt: SEQ ID NO: 20

caaaaataacagcaagaaaagcggaaagaccatcgcaaggtggaaaggat	50
tataatggcacagcaaagtcgcacagagcactacagtatagcatagagtg	100
ctaattgagttgataggcccaatTTTgattatgccttctTTTtcatacacg	150
acgccagaggacattattacattacagtagttcgccgctagatgacaaac	200
gacatccttaccgatatgagatgtgcaaagctacataatggcaacaagcg	250
ttatgaacagccttgtctTTTtacgaccacagaaaagccgtattagagctct	300
tcagctgcaaaatTTTcttctaataatgatgcaaagccatcaaaaatcatg	350
catagttatgaaatacctgatgaaacgcttcgagttcgtgctcaagaaat	400
tactgaaaggttaccgagaagaaaaatatctatgagacacgataaggccc	450
cttctgaatccattgtcctgggcttgttcattctatTTTaccacttaaaat	500
tgatcctTTTcaaaggaatTTTTTctatTTTccaatagtatatTTTgtacaa	550
aaactacaaaaatggataaaaaataacagtaatttTgtgactactgtaaat	600
atcactgatTTTggattTTTgtaatgagtactgctcatgcccattgccgatgc	650
aagtggatcataaattTTTactaaacgatattcgataatgcgccaagcctt	700
tataagggaactcaaaataaccatattggacagtttcagaaggccaaataa	750
cgatcaaggacattcactcatgtTTTTTcaaaggcgaagagtgtaaaattt	800
tcttctatatagttcgaatTTTTTatcttataaatttcagtcgtcatttt	850
ccacattcgaactcaaataatgataaagaacgctgcagtaatggcttaaa	900

FIG. 25

26/51

## Utility3

nt: SEQ ID NO: 21

gcaagtatggggtagcaagctgcttaaagcttcttatcacactgtaacct	50
ctgttaaagcagctgcgttggtcttgcaactctgcacgtaaatacgcacgg	100
ccaatgtaaaaatagcagtcactgcaggcctcttggtattgtacccgtaat	150
tacagcaattgccattagattttacgaagtcatttaattaaatggtttgat	200
ttatctccttaaatgcttccagaaatgggactagacttttttctactcaaa	250
cctgttcacaattattgctctttttcaattataaggtaaacaaggccatc	300
tatcagcaacacagtgctcgcatTTTTtaattaaactatataaaaccaac	350
tatttggtggttgcgacttcactttttgttgaattactaccaatcattaa	400
tattgaagatgtgagatcatagatttattggctttgggcatctcaaattcc	450
caagaggtcattttaaccaacaacattttaaaaagtagatttgtctgcct	500
cagctatgagatgcgcatgtccctagcatctcatatctggttatatatt	550
ttttccacttggtgaatgttgaaaaaacaccactcgtccaatttatcag	600
tttgcaggtctaattgtccttcctgttatttaaactgtatattgtaagca	650
tgtcttatcgaaacaacttactcagttgtccgaaaacaaaactgcaaatt	700
ctgtgtgtattcacgtactagaatcctgtcaaattggatcttgatttaag	750
cttttatagcaacgaactttgcatactaagttttttttgttaaccggaac	800
tgccaagaagcattcagtaaaatacatcttcatcatttactgataatact	850
cattcagactcatatcatactatttcgaattcattatacatcctcaaaaa	900
ccatattcttcagttgtaataaaagatagagcctgcatttgattcgattt	950

FIG. 26

27/51

Negative1

nt: SEQ ID NO: 22

gatttaatacagtagctttcttcgctaggatctatatgcgaatatatcac	50
atatgtaaattataagctcatcgcaaaaccccccccccaattttcaa	100
taatttttccactaatcttcaaaaacaaatggggtaacccgtacaagagtt	150
attaaaacccaaaatgacaaaatcgcgacaattcaatcctacttaattag	200
caataacatactagcggtagagctactatcacatggtgaaccttgaatgc	250
tcaattcattgtactcaatactgctatcaaaagaaaaaaatgtattaat	300
tatattcttgtcaaaatcaattttacactataagaggaaaatgttcttca	350
gtcctagtaacatttagttttctcccttgctagagactttacataatatc	400
ctagaaggtaaaattcgataatacagcagtaaaagtcgtatattggtagca	450
atccttggtgacgctgacttttttttttgaattttattgttttagttca	500
tgataaaaaacttcaaataccttttaatactggtagacagagaaaacaaat	550
cgaaacgaaaatagagaactacgaataaaaaaatataagtggagaagatc	600
gtcactacgcattaaacaatattgatcgctcaatgccagtagtgcgcgta	650
aaagtttagtaacttaacgatttaggcacaaatttgagaaaaatttcgccc	700
tgcagtaagtatgttattcagtagcatataaagctgaggttttatgctgg	750
caacgttcagatttttttaggttatcagcaatgttaaaatattaaatagga	800
tacttttattgtttgagaccacctcaatgccagatatgttaaacgcttt	850
tttctggagtgaggatcatagaaaaaggctcgagtacatcaagcactta	900
aaggttcaacactctactgttacttctttaagctaagctattcatacata	950
atagtcacatcaaagtgg	967

FIG. 27

28/51

## Negative2

nt: SEQ ID NO: 23

gcaatttgcagttcaacttttcaatgatgatttagaatgatctaactg	50
gaagttgaagtttttcaaaaatttgctgtaaaatttgaccgatttgtgag	100
attcttcctggctgtcagaatatggggccgtagtatattgtcagacctgt	150
tcctttaagaggtgatggtgataggcgttgagtatgtgtagtgtttgacc	200
cgaggggtatggttttcacaagtactgcgcactgtattgtgaaagcagctt	250
cgggggtgcgtgattaaaaaatgcgaccaagaataaacagggtacatcataa	300
caagggccatttgaattgcatttatcaggatttgtaaccttggttctaaag	350
aggcatcgatatagtttaagttcattttccacccaatttgatgacggtgtg	400
gaccttaacctattgtcttgaaatttaggttatctcttagatatcacatg	450
tgattaccccagtgaaacgcgtataagcttacagaaaggaaaaccggttg	500
ctcagtcaaaactggttgagatttgggctcccctgaatatttgagacatc	550
cctaaaatgaagagatatatacagctaattttgaatgaaaattttaaatt	600
cgcaatgaacagtactagagatgagcttttgaagtcctttcaaattat	650
gttcttccagttgatattttttatatttatataccagttaccaaataatac	700
ttgccatacatttacctttttgaggttggttcaacggaaatccagtgatt	750
tacacattcttggaacccatcgcttataatacgaactaatttat	800
aacaaaggcttttggaaggtatccctactttttacgacgctaaatcatga	850
tacgaaacttttaggaagattaacagtcactccataaaatcagaaagtatt	900
cgctaatagtggaaagaaatgggttatataaagatggaaatatcttgaaag	950
agacagtttaaccgaagttctgtcaaagt	981

FIG. 28



29/51

DR1s

nt: SEQ ID NO: 24

ttcagaaaagcaaggaaacagtactatcgtttagaatgtagaatgatagg	50
ttgcttgctaattctattatggcacgaatgatacacccatattttcaaca	100
aaatcaatacccactagcatcattgagccaactatttgtcaatgcaacca	150
ttaccggtacttcatcctgatttaacgagtctacttttttatcacgtcaa	200
aatttacttgttttcctgtaaacccgaaataaaggcaaaaaagacctggg	250
tgcaattacgaataaatgtacaataatcatcctgtttgcatagtaaactt	300
ccagttagagtcacacaacgcaatgaattttgacagttttctgtgcgata	350
ttctttggtaaacgtaaagaacaggcaacttttggtacaatggattctag	400
cccatatgggttcatttctgggtgcattcgcaaagtcagtatttgtctagct	450
gtgttttctggctgagagacattatgatgttattcattgttatggatatac	500
tctgtagctcatgctgcttatttctccctaaaaaagtttttctctcgaa	550
tacattcttgaccatttcatagtgaattcttgtacttatttaaaaccaa	600
aaatggaagtattcatatcatccccctatcaaaaacactcaataagtttcg	650
aattattcgttcgtctaaacagtggtccaatactcaaaggggtattcaaga	700
cggcacaaaatcagcatcttcccttatccgtgttccagaaataccacgct	750
aaggtttttccctcctacaatccataaaatcattaaggaggcagcttgaaa	800
aatcttg	807

FIG. 29

30/51

DR2s

nt: SEQ ID NO: 25

tttctttttccctatctcactgggttgacagaaatcagtgtgctatcat	50
cctaccatatgcgctaaacttattgtctttctcctcctagagatgctgta	100
ttccatgcatattctgaacgatgggttggtgtttttatcaagcaaggta	150
atcacatggcgtggcttgctccacacatcagtagaaaacgcataccgcag	200
cggaatccttaaataataagtgattttactgttcatcaactacaatcgga	250
ctctttcacaaattacccttcttggtttccacatttactgttaaatagaagg	300
gatgtacagaaggcttaggaaaacctgtgctgaatactggatggacactg	350
cattcccacagtgaacttttatagatacactgtcagttattttcgaact	400
ttcatcaagttgctgagtttttagtatccctttgccttagctatatgtttg	450
aatgagcaaaatatttgcaatgtctctagctttcttgaaatattgggtta	500
tattgagggccttggttaagatttcaaatttcactttgaaatactcaggaga	550
aaaatcatgctcttttgataatttggtgactaaacatacataaaacagtt	600
taattttgggtggtaatggctgtgtgactagctatagaaagaaaaaaatt	650
aaaaaaaaaaaaaaaaatcaagtagttcctgcactgcgacgtccattata	700
gcattatgaattgggtccctgatttacgcatgcgataaactattttttagcg	750
cagccgcatatt	762

FIG. 30

31/51

DR3s

nt: SEQ ID NO: 26

aaaatctcaaaattcccaatattcacatagtctaaagtaccgatagcaac	50
caacatatataaacagtagtattttacgaagctgaattgcaagattagt	100
agaggagaataaccggataatTTTTTTGGATTACGTTATTGTTAAAGGCT	150
ataatattaggtgaaacagaatgtcctagaagtttttttctttcatgtta	200
aatttattgattcttgcgcttcagcttttataaaacataagaactgtttc	250
ttcacgttaacttcttgtgccacatatataatgatgtactagtaatatgggt	300
actatttggcagatgatatttgaatttttattcaagacggttactgtttct	350
acgattgatattttcattcctggatatcatcttgccagatcacttacaat	400
ttaggccgcgcctgaattgaagagtacttcaatacgtagtgtactgtcca	450
aactctcttccaaatTTTtaatatTTtagctgggggttgggtaacaagtgag	500
caagggaaaaagtgaacattTTtaagaagaacaataaaatagcaagagatg	550
gaatggtaatgcttggctctcgagaagagtagcataaaacgagacttggt	600
taaaacaggatatgacatacttcaattcagctttccctatcagccgctcg	650
agcagttatataggtgtgttgccggagtaatttggcggaggccaacagt	700
gctaggcggcaacgcctggaacacgcgcttaaaagttctggaagggttcgc	750
gaattgagaactgctcaggggcgaatacaggggcggccttggcggcaggg	800
gggaggcctctgtgaagttagttatataagacttgctgtcatcgTTTTT	850
tgatcccggcaggaactatctttt	874

FIG. 31

32/51

DR4s

nt: SEQ ID NO: 27

aatagatgaattatgtgccgctggacgtttatagatagcataagcacaat	50
gacttaaagggttataataactcattgatatcactctgattataaaatcgta	100
atatgcgaataggtgaactaatcggaataaccatacgacacttcaagct	150
tcaattctatttcaactgtagtgctgctagtgagaatacaaaagtagc	200
atacgtgatgtgcaaaaaatgcgctacttatcacacaagtaccttgcgca	250
agaagggtactctaaaccggggccatcgccattaccagacggagatgtatt	300
ctttatgaagcaataattggaggtgtatcaagttcgaaactgctgatgct	350
atggatttacatctttcttatgcacaaggcttgcttggtttctgagtag	400
ttagtttttagatttttgtcaagtctggggaagttaattcgagcaaaat	450
taacggcacgttattctaatgcatatgttggtcatatattcttttacaaa	500
gaggtttggaatgatgtcaccgatgttagaatgttaggagaatttcatgt	550
gaattttagtccaagtgttgaagttctctctgcagttagggcacgtaca	600
tggcaacgatatcgtttttgatgtattaatcttagtaggcgttgagtttg	650
tatgttacttttctcaggtgatgaagcgtgatgacgatgacaaaaatggg	700
ttataatagggcgactatcatcatgcgtgattgatatttaaccaatgtc	750
ttgagtacatcaactccagaaaatgggtcatttatatgcctagcatgt	797

FIG. 32

33/51

DR5s

nt: SEQ ID NO: 28

tggcaatgatactcgttattcgtaatatcagtcctgcaaggtgctgtgat	50
ttctctatTTTTatattgcctattatTTTTTcaaagatttgagccgTTTT	100
aaattgagtatgcaatgagtcTTTTgaatcaaccgtaaggcagttccata	150
accactgccacgaatacgtttcactaccttgaagaatctcctaattgtaggc	200
cgtattcttctgcacttagttctgacgatgtagacatctcattatataaga	250
gcataagcgccctgtttctagaatcatttcttctgacccagctTTTTtgag	300
ttatttctcgcggtattttgaaacatttctctgagcttgacgtgaacatcctt	350
atatttcatgacaaactcgatcattggaacatccctgcctcgatttttaga	400
gctagtatcaaatttcaatctctttgtgatggagccccgctcctatttca	450
aaagagaagtttcttgtatgcatatgttattgaagtctgattatagcaag	500
tgcaatgtcgtctcaattattttaactatTTTTtagccatacatgttagtt	550
atcctcaaagagagcctccagactgggaagcagtgtttgtcatttcaa	600
aagtagatttcacagtttgtatgattttcgaagccaggattcattgggct	650
ttgagtaaagagaagccgcgtattacgaacagcttacgatattgtaaaat	700
attcccttattgtggtgccccaatggatacatgccagagaaatgtctgtg	750
aaattgaacaattacaatgacgagagcaagtaatccggcggccttgtctc	800
tctttcac	808

FIG. 33

34/51

DR6s

nt: SEQ ID NO: 29

agataccgtccttggatagagcgctggagatagctggctctcaatctggtg	50
gagtaccatgggacaccagtgatgactctagtgacttgatcagcgggaat	100
accagtcaacatagtggtgaaatcacgtagttgaaaacagcttcagcaa	150
tttcaactgggtaagtcttcagttggatgagcagcttggaacatatagtat	200
tcagccaaatgagctctgatatctgagacgtagacacctaattcgaccag	250
gttaactctttcgtcagagggagataaagtagtggtggctggggcagcag	300
cgacaccagcagcaatagcagcgacaccagcaacaattgaagttagtttg	350
accatcttttttcgattgaactttttagatcttttttagtgaagatgtgag	400
ctcactcgaatgtaaataacaatgccaaattgtcggaaagagttaatcaa	450
agctgctctatatttatatgccgttttttaataagcgacggacgaacagata	500
aattggtgaatagctatcttactgctgatatttctcttacttgggctccc	550
ctatcccatactcttcaccactacaaatatgcagttgccctttcttcaac	600
aatgcttttttttatagatctcgtatacggatccgcgcctttgtactacct	650
atatcttattatgatataacaggagcacaggaatgttcggtacagggat	700
gataccttt	709

FIG. 34

35/51

DR7s

nt: SEQ ID NO: 30

ttgggacggtttttgcactaagaacagacgagtttacggttatcctcaac	50
aagcaagcaagtatttgctaactagatgccattccgaatcattactcat	100
acgttactattgagagatgttttacaatagatgagaagaatacaatgtcc	150
agagctcctggtagtagagtgcatattccaggtcttattcgaatcata	200
tcataccgtccatttcaacaatggtgaaatgtggtccacatatatcagaa	250
atcttaacatttagtgaggagagccagtagaaaaatgtgcgcaagcggaa	300
agaagtcatcacagacacgtttaacaaaacaccaccacagcagctttgt	350
ctcttgattctgatcagtttgccatcgaagaagcaaaattgtggtgttat	400
ttttttcaaacaaaacttttttggcaacagcagttttcttctggatattt	450
gtactttatcatccaaccgatgaaagctggtttcctgtcaacctacattt	500
aaatggcccgtacttcttcaaaaccgctagataagcaaattaacccaact	550
tttgagcgtcctaaattccccttggctcagaagactcgtaatatgggaa	600
gtttaagtcctaccatataatcaaattggaagctttctgtgttcgaatgg	650
ctattctaaccgctgggctattaatcagaggggaagtgaaatgaccgaga	700
cgtattatacgtcatgttgacatcaacaatttaaggaaaaaaataaaaaa	750
aagcaatgaaaaagggtttttttaagttgaagacccttttcaaatatatg	800
ttgctttgaattgtatctaccgtctcggtttcttctgctttaccgtttttt	850
tttgcccttctttagatatgtcttttatgcttgaaaggtccggc	893

FIG. 35

36/51

DR8s

nt: SEQ ID NO: 31

ttcacgcctagaaaacagcggttgctgagaaaaaaataaatcatcgag	50
aagaagtatgtcatcataggatgttccattgtaaggatgatgtgaacat	100
actcgaacaaagaatgtatagagctgaatatttctcctttaaatttcaa	150
gaaaatgagaaggaaaatctcaaacagaaacttcggttctttttctcaagt	200
aagcaaaagcttattgagacaaagcggaataactacgatattaataacgt	250
tgatgaagctcgaacaaagttagcgctcggttatgcttgccatataaaga	300
tatatttgccttacattttcgttgaacgtagaatgatttttgctttta	350
aaattttttggttgttctttcagtgcttcttcaactttgatacgaaagcaa	400
gtgcattagtagacaacaagaactggccacaactatactatactcattttt	450
cttgcccggtgttttaaagtgtttcatccacagcatttgatgggatgattg	500
gaagtgagacggttcgagaaaatccatattttgagtcaagaattcagataa	550
tatactgagatgattaggtatggctgggttctacaaaaacacaaatatcc	600
ggctagcaatgatcactgagcaaattaaagcggttaactcactcattattg	650
tagcttatgcgtttctcctcctctcttttttctcgaaccggagtggaa	700
gatccaataacgtaataattactgatgttggtatttaaagctggcaaaaata	750
acatgaggcgtaaaaccgcactgcggttaagatgaggg	788

FIG. 36



37/51

DR9s

nt: SEQ ID NO: 32

tgaaaaagaaatatttgcgaccttttttcggttacattgatcgtgaaattt	50
taatcaaagataatataaggacgtgagatatttatctttttacttgaaat	100
taacaatagaattgcgctaagcggaataagagcttttcgtaaaccctttcta	150
tttgcaccattgcgtcaacgtataaaaatggtatgacctttacacaaacgc	200
atgcttataatcttatgtttttcataggggtgtaatttggttgatgacgta	250
gtctaaatttgatgctatctgcaattgaggtacatataagaggtcaattt	300
cgggaccaacccttttaatcgaaaaaacgtaattcactagggcaaggga	350
gaacttagcagctaataatcgtaaaccctttcataactaaaaaatgcactta	400
ccatcaacaaaaaactcaggaccaatttccaagcttttctaggtgattgc	450
ctataacacaaaaagattcgctcatatgagatttttacatgtaatagc	500
aatttgttccgatcagttgaagggtcatcaacgcacggcaggtacatccac	550
acctatcacaaagcccttcaataattcacctacgtaaagttataccgaaa	600
catgcaaaatccatgaaaaattctgtatgataacgatcatatccttttgt	650
attggtggtacgatgctcaaagatagttattggttgcacctgaggcaaaag	700
cggaaatgaaaaatccagatggggccaaaagcagaagtattgtgtacaac	750
aattgcttcagcagtttaccaaaccgtttccagcaatcatcaaaagttg	800
ctttagccacatttccgcaagatat	825

FIG. 37

38/51

DR10s

nt: SEQ ID NO: 33

attgaagctgggtgtggaagattttatttgaagaaactaaaacgtaccctg	50
tcatttcctgagtcctcctttcaacttagtgtgaaagccgaacaattataa	100
tcctcggtagacaacagattttattgtactaaagttactcttcctgttatac	150
ttccttgattttactgttatagcaatgacccaccgcaatcaggagagccg	200
ccgtatggaatagcataccaagtcataaaatcgtcaacctattaacgggg	250
ttcaggttctttttcagcgtagtagccctttaacaagcgtgacaaagtt	300
gacactcagagaaaattcaggattttattgtaatccagctactcatcctta	350
gatccgcttgaggcatgggtttttttcaccttgagaggctattttgggta	400
agccaggaaggctgaaaaatcccaaaggacacagtaataagaaattggt	450
gttggtgtatgatgcatttagaactcaaaagacgagtttctgaaaatgct	500
tacaatactccataggtaacatgatttttttattaaaaaagtatactgtt	550
cctttgggtaaaaattatgcaacccttgagtggtccgatgaagataagact	600
acgaaacaatttgcggtaaattttttctgctatttgacatttacacatgct	650
ccaatccattaccctttccattctcgtaataaaaacctcgaactgttattt	700
catatttacatctagacgggtatcggcctcaacaactcaaacaaaagta	750
aatagaaaagagccagacctatcgc	775

FIG. 38

39/51

RC1s

nt: SEQ ID NO: 34

gctgcaagtcaatctacgggaaagaagaaatTTTTTaaacctaataatgcaaa	50
ataagcttttcttgaaaataagatTTTtcggcaataaaaaggtaaatagcag	100
ccaaaaatcaaaatacttcagaagaagtcgtagcgaggactgctaagggg	150
aagcggatttgaagatcctttccagaacaagaaggagccgaaagctgtca	200
ggaactgttcctgattTTTTtaggaaaacaattaataggtatctcgtctag	250
agtagtatctcgagcttccagaagttgcagataatcaaaatcattgtttt	300
atccctTTTTtagattacagcttagaagagtagagagcaagtttactga	350
aacggttccttgtttacaataatattcctaacaaactttacgaattagga	400
tgcagcatgattTTTTtatattgcttcacttcctaaagtatgaatttttat	450
ccgtagtcgcaaacaaaacagctactggaaatctgcagcttggtaaaaac	500
cggtagttccgaatactcctcgctccttgagttgtataccggttaaacttc	550
ctaggggtgcatgtgtctggcccaattggcccacaaaatctggtcctatt	600
gacggtttcttttgattttcagcatcttcctctaagaaggacagaaaat	650
tatgtaatatatgggagaaacggcctcccaactgctaagtgtccccggca	700
gcacgagtaagcaaaattcaggcaaaactattgcattaagaagccgtacat	750
aattcagcgtgatatgatgaaatTTTgttaattgcaaatttttagtacgat	800
ttggttgttagtgtgtgtttatgcaagtaattattgaaccctaagtagtt	850
actgtcttcttttgctgtaattcgtggattcacg	884

FIG. 39

40/51

RC2s

nt: SEQ ID NO: 35

gtccccgtttctcatttttgagacatgatctgaacaaggctgaaaacagc	50
aatctttttcgataacttttgcaaaaatttcaaacattggttggttgatg	100
cagccaatttttatagggtagagagcttaatgctttacatgtgctttatt	150
ttcgggtactttccttaaagtgtctacattatctctcaggacttgaatgtc	200
ttcgggtgaattactataaaaatcttgagttttctctgaagtttaataccta	250
agacaatagtggtgagtgatgtagttcacgtgtgtgccactggtaataat	300
agagataactatctcagttaagtttgaaaaggtaaaaaatagtttaagta	350
gtcattttttgcgacgggtcattcttctctgatgcagttcttttagactac	400
ctataaacaccattcttacggaattataatggaaataaaacatcagtacg	450
tgttgctgtcgggtgatagaggggtaacagaaaccttaattgaaaaattagc	500
acagtgcataatttattaacatgattgttttctgtgggaaataagaaatt	550
tcagcaccagtaaaagacgagaaatatagggcacataaatgcgctcttac	600
tcgtatgttccaggatgaaaatgttttagggcatcaagtattgccgaaagg	650
gcaatatgctttaacaccagaaaatccactgtatactcgttacgggtaaa	700
caaagcaaaacgcagtgcggtgataatgtttctaaaatctctgcacactgt	750
tgaaatgcggctctgatacttttagcc	776

FIG. 40

41/51

RC3s

nt: SEQ ID NO: 36

ccagattgcttacaaaagaatagcgagccaacatttgctctgcctcaggc	50
ctcttggtgctgcttgaagactcatcttatatggcttttgatgtcatga	100
tttggttcttgctacattatgtgttgatattaacaaattgattttttttt	150
tttgcgatagcaagcagataatgaaagagacaaggacttggaacatccga	200
taagactgcgcgatatcgatcttacagtccttcccttggtgcatgactt	250
tcggaaaagcatcctcgtcgactggtagtttgctgtctgtcacgtgctga	300
agggctctgatacatTTTTTTaaagataagagacgggggttacccttcgga	350
ggactaagcgagatctccaagtaaagatctcgcttatcaagaaagcagcc	400
aagtgtggaacgtcctTTTTTTTgggttcaaaaagatattcaacagttta	450
cactgcagctttaattgcctcaaaaggatatcatgaggtgatctagggtc	500
agaagggaaagattacagcatcttgagttgaatcacatctgcaaaagggtg	550
gtattattgacgttgctcttcccttaatggaaactcatgggggttggaag	600
gaggtgcggtaatctatTTTTTTcgaacacaaaacctaaccttgaaaaga	650
aactgtccaatttcattgaacttacctcagaacggggccggagtctttgct	700
ttcagttctaactg	714

FIG. 41

42/51

RC4s

nt: SEQ ID NO: 37

ttcgcgtattccttacatcttcgaagagaacttctggtgtaagtataataa	50
atattatagctctatcgaatggtgcaattatctaccaaattctcaatagg	100
aatccataatactacatacagataactaattcttagtatttttatacttat	150
tatttctttttttattacaccagcaatcggttgcaaattatcttctgataga	200
atctctgaggggtatcctaaacttatgccattttcttggtgactgtaaactcat	250
acttggtatggtgtgcatttagtcaataatcggttcttggtccaacgattac	300
atgtaaatagaaggagaaataattatggtaaactcatgcggcggtcctttt	350
ggtgatgcagtatccatagtcactacataacaactcttagtcaccttgat	400
tgattcaccacataatcctgcagagcccgctatgtccttaactctgcgcga	450
taactctcctacccctgaattttgagagcgccatagcaaaccgataaagc	500
tggcacaattaaaggtatcggtgttgtcagaattaggtgcctcctgcttt	550
tttttttttctgctcttatatccggttatatccgaatgatttttatcgct	600
tgtttaaaaaatactttcccgatatatatatatagtcctccctttaaattt	650
gtttccggttaagtttttaacaccaataaatgaaaagaaatgactacggtg	700
atgaatatgagccgcgcattgaatcaggttatgtaagtatcagaaccct	750
aattatg	757

FIG. 42

43/51

RC5s

nt: SEQ ID NO: 38

ctcaagaacggtggttgggtgcatcaaaagttttcgactgcttatttggtc	50
ggaaatataaaaaactcgatcctcttatctaagcagtatacattcttcttt	100
ttgaaatgaatgtactccgtaatatcttcttatttggcattttcatcctt	150
aacttttgcattggctctgaactagtcagatagttgcccttttcagcaaac	200
ctcttattattgaaagcatgggtgtacatccgttatactattatattataa	250
gaaattgggatgccaatTTTTTtgcttttggtttgcctgttttccttctt	300
ttcgcaaaagtaattgcagatttaatagcaggatattataaccgttggtaa	350
aacttaaggattttatgaacaatagcttcaagtacagcattcatagaacc	400
aactactaaggatgaaactagtatgtttttgtcaaaatattttcttgacc	450
ttgctgtaacatcaagatctgtttctctaagatattaaagttgagtaaaa	500
acaaagctgatatgagaaaaatacgtaattgctccacataatacgtgggt	550
cagacataaaggtagaataacttgatacagaagagattattcggtactctt	600
gatggcgtgcttgaactgggtgcctcttaacaaccggtaatatagtcagat	650
gagtcactacgagtgtgtgtagtagcaagtgttttacctacgtggcagta	700
agagtagctctatggttgtgtaatagtggtgcttattcctaatactctga	750
agtctgaagcggtagcagttgggtctgtatcatgggtcaaaggagcaa	800
acatatcttctgaagtgaccgcaaatagtactatgatgtggttggcaata	850
taacttaaaaggaaataaccacaaggaattgcacccatgta	891

FIG. 43

44/51

RC6s

nt: SEQ ID NO: 39

tttttatatgggggtctggcgcttcgggaaaagagaggaaaacttgtaac	50
tcaatatatctcgatacaacattacgttttgtaaatttatcacaaaagcc	100
aatgatgatatctctcttgcaagttatcgaaacattgattggtaatttgt	150
ttgaaaattgttaattttattgaatatttcttttgcaaaagaaatagtctc	200
agcgaaagctgggttacaaaatttacatcatgagtttacgggatttgtaaa	250
tacgctttttgcataaaaatactttgccgtttcccacccttgcatattca	300
cttactcccccttcatatactctatgtaatgatgattaagctttggccg	350
ctaagtctctcaattagtgttgatttttggtttttattcatatgattcttct	400
ttagtgaagtattgatcaattacgtgagtcagctttttgaaaacccatt	450
tggaaggaattaggaaattattttgcttactacgaccactaatttaccgc	500
catttctgggcctttttattgactattttgaccatgtgctcgactagaag	550
aacggcatcataatctgctggtagagtttagtctataatgattgttgaaaa	600
taaaggcataagagatattccacctaataattcaagttattgactttatta	650
tcaggatcttagtatccttttttggttaagtcatttcaatgaactaggtc	700
tcgcaaactttttgttcgaaaagcggtagtgcatagttatgctaactctg	750
gatatatggcataaacggtacaacactagcccatttttttgggaagtagtg	800
agggcagctagactgtatgatgaatattcgccctgcatactgagttttt	848

FIG. 44



45/51

RC7s

nt: SEQ ID NO: 40

tgaagacaggattcaaaaccgattaatagtagcagaaactaaaaaagtac	50
gaatattagtaaaattcatgttcttgaatcgagctactatctttgtcggg	100
agggtaaacgattataactcaaaatgactggaactggtgattattaattt	150
ttacgtttcctgtgccataagcgggaagataagaggatagaagaaaagaa	200
aggcggcacttggcgaactacaatggcgattatattcatggcgattatat	250
tcatacaaaggtaatggaggcctcggataatggacaatattgagaaaatc	300
cttatgcttacttctcttaataaaaaatagacacagccatttattatgcg	350
taaaaaagattaccacttgtcttcgatgcgtgctgctgccaatcaacct	400
tttgagcggaaacttcgagctcgcaatgcgtctggaatggtgctagagaca	450
gtcttggttatctgtgacatgtgtttcgttcaggcgtgtgagcatcttct	500
tgttcgatttcaaaattaccgccttgactcgtgaaactggataattcggt	550
ggcgttttcatataagtcgtctgatggcgaaaacttttcctttacttagc	600
atacagcaaatatccccatttgacggatttttgaaaaatgagcccgctaa	650
cccagaatgaactgcattaccaagcatttatgtaaacgttccgccaccat	700
ctttggtaagggtatactattatgttctggatttaagggttgattcacaatt	750
tttcatcaccaaaatctggtggcatgcctagttgtctggtttcaggcaat	800
ttagcc	806

FIG. 45

46/51

RC8s

nt: SEQ ID NO: 41

aatttgagagcggatcttgcattttttatttatcatgcttatgttttttc	50
tttgatgtaagaagaagcaagtaagatatgtgaatatcttatcactaatt	100
caaataactaagagagctcacaacgacaatttgtgacagcatgcgaagca	150
aagagcagtgataccagtatctttcatccagtaataacatacgactgatg	200
ttatagttaaatgttacattttgagagacttcaacctctcgaaaccaaga	250
ggttggttttaactctggtgacttcaagaaggggtgggtaccttttataaa	300
gcttgagacgaagcaatagtcagtcctctgtataacaaggagaccacctca	350
ttttccagtaactcttgaggcatgtcggatggtttgccttgaataaacgg	400
cagtcattataatgaatggcctgtactttcaaaacagtcctggaaacagaa	450
atccattgctgaggtaccttttagtagcactttcgttagtgaagggttaa	500
ggtagttcttatttactgcacaagagtttacatttaaccactctaatag	550
taactgtagagtggtttaactgttaggtgatctgttcattccatttttc	600
gtggtgtatctcaagatgagatagcttagcggttgctacatacataaatct	650
aaacatataaacacctgtgtaactcgtaaacgtctgggcttccatgcttc	700
taccatttagaatgatgtagaccatttattccaagaggataagcaccctc	750
tgtgattcaaaat	763

FIG. 46

47/51

## Utility1s

nt: SEQ ID NO: 42

ttaatggacttttttagtgtcatcgaattttatgtaatatataagaaggta	50
gaataatttggcaggataatgtgttagcaaaggaggaaatcgaatacctt	100
taaaagagaaaaaatttttttagctgcttaaatctctgtgttataaccaccc	150
gatagattttgagttatgctttctaattgatctgactgcgaacgttttct	200
ttatgccatctgaattgtcaggaacaaagaagaaaaagaaaagtttttaa	250
aaaatctgtggtcgtgtgtgatgtacctttcctttacatgcattaatgcg	300
ctctgaaatgtggtacgatatccttacagagaatatattttctgtatctc	350
gtgcaatgttgaataacctatgaaggaaagtacccatcgctcaaggtaag	400
cattccaggagggtcgccagaaacttaaaactagtttttagcgacagatccg	450
aaaattgatagagacattgaaaaaatcactactccgtccttttttagtgct	500
ttctcaatgcataatttttggtgcacgactaaaaaattctagaacactata	550
gttgcatttttttgggccggaagaagaaaaaacgcatgtaactttaatgtca	600
aataaagttttcacctagtaagcgcgatacaaaaaaacacagaaatagc	650
cataggaaagtgaattttgtcagccgactaaaattaagggttagcttaca	700
agcagcaaaaaatttgacatcgacgggtattccctgaaaaaggagcaggc	750
aggtgctgtatatatttttttcggttcctgcctcttacatggcgctcggtgta	800
tcttaaataactaaagtgagctgactacccttttgagtgccctatgtgacc	850
tctgatctcgaaagtaacaaga	873

FIG. 47

48/51

## Utility2s

nt: SEQ ID NO: 43

ttataatggcacagcaaagtcgcacagagcactacagtatagcatagagt	50
gctaataagagttgataggcccaattttgattatgccttctttttcatacac	100
gacgccagaggacattattacattacagtagttcgccgctagatgacaaa	150
cgacatccttaccgatatgagatgtgcaaagctacataatggcaacaagc	200
gttatgaacagccttgtctttacgaccacagaaaagccgtattagagctc	250
ttcagctgcaaaaattttcttctaataatgatgcaaagccatcaaaaatcat	300
gcatagttatgaaatacctgatgaaacgcttcgagttcgtgctcaagaaa	350
ttactgaaagggttaccgagaagaaaaatatctatgagacacgataaggcc	400
ccttctgaatccattgtcctgggcttgttcattctatttaccacttaaaa	450
ttgatcctttcaaaggaatttttttctatttccaatagtatatatttgtaca	500
aaaactacaaaaatggataaaaaataacagtaatttgtgactactgtaaa	550
tatcactgatttggatttttgtaatgagtactgctcatgcccattgccgatg	600
caagtggatcataaattttactaaacgatattcgataatgcgccaagcct	650
ttataaggaactcaaaataacccatatggacagtttcagaaggccaata	700
acgatcaaggacattcactcatgtttttcaaaggcgaagagtgtaaaatt	750
ttcttctatatagttcgaatattttatcttataaatttcagtcgtcattt	800

FIG. 48

49/51

## Utility3s

nt: SEQ ID NO: 44

tctgttaaagcagctgcgttggttcttgactctgcacgtaaatagacgacg	50
gccaatgtaaaaatagcagtcactgcaggcctcttgattgtaccgtaa	100
ttacagcaattgccattagatttacgaagtcatttaattaaatggtttga	150
tttatctccttaaatagttccagaaatgggactagacttttttactcaa	200
acctgttcacaattattgctctttttcaattataaggtaaacaaggccat	250
ctatcagcaacacagtgctcgcattttttaattaaactatataaaaccaa	300
ctatttgtgggttgcgacttcactttttgttgaattactaccaatcatta	350
atattgaagatgtgagatcatagatttattggctttgggcatctcaaatac	400
ccaagaggtcattttaaccaacaacattttaaaaagtagatttgtctgcc	450
tcagctatgagatgcgcattgtccctagcatctcatatctgggttatattat	500
tttttccacttggtgaatgttgaaaaaacaccactcgtccaatttatca	550
gtttgcaggtctaattgtccttcctgttattttaaaactgtatattgtaagc	600
atgtcttatcgaaacaacttactcagttgtccgaaaacaaaactgcaaat	650
tctgtgtgtattcacgtactagaatcctgtcaaattggatcttgatttaa	700
gcttttatagcaacgaactttgcatactaagttttttttgttaaccggaa	750
ctgccaagaagcattcagtaaaatacatcttcatcatttactgataatac	800
tcattcagactcatatcatactatttcgaattcattatacatcctcaaaa	850

FIG. 49

50/51

## Negativels

nt: SEQ ID NO: 45

gctaggatctatatgCGAATATATCACATATGTAAATTATAAGCTCATCG	50
CAAAACCAAAAAAAAAAATTTTCAATAATTTTTCACTAATCTTCAAAA	100
ACAAATGGGGTAACCCGTACAAGAGTTATTAAAACCCAAATGACAAAAT	150
CGCGACAATTCAATCCTACTTAATTAGCAATAACATACTAGCGGTAGAGC	200
TACTATCACATGTTGAACCTTGAATGCTCAATTCATTGTACTCAATACTG	250
CTATCAAAAGAAAAAAAAATGTATTAATTATATTCTTGTCAAAATCAATTT	300
TACACTATAAGAGGAAAATGTTCTTCAGTCCTAGTAACATTAGTTTTCTC	350
CCTTTGCTAGAGACTTTACATAATATCCTAGAAGGTAAAATTCGATAATA	400
CAGCAGTAAAGTCGTATATTGGTAGCAATCCTTGGTGACGCTGACTTTTT	450
TTTTTTGTAATTTTATTGTTTAGTTTCATGATAAAAAACTTCAAATCACTT	500
TTAATCTGGTAGACAGAGAAAACAAATCGAAACGAAAATAGAGAACTACG	550
AATAAAAAATATAAGTGGAGAAGATCGTCCTACGCATTAAACAATATT	600
GATCGCTCAATGCCAGTACTGCGCGTAAAAGTTTAGTAACTTAACGATTT	650
AGGCACAATTTGAGAAAAATTTGCGCCTGCAGTAAGTATGTTATTCAGTA	700
CGATATAAAGCTGAGGTTTTATGCT	725

FIG. 50

51/51

Negative2s

nt: SEQ ID NO: 46

ggaagttgaagtttttcaaaaatttgctgtaaaatttgaccgatttgtga	50
gattcttcctggctgtcagaatatggggccgtagtatattgtcagacctg	100
ttcctttaagaggtgatggtgataggcgtagtatgtgtagtggttgac	150
ccgaggggatggttttcacaagtactgcgcactgtattgtgaaagcagct	200
tcgggggtgctgattaaaaaatgcgaccaagaataaacagggtacatcata	250
acaaggggccatttgaattgcatttatcaggatttgtaaccttggttctaaa	300
gaggcatcgtatagtttaagttcattttccaccaatttgatgacgggtgt	350
ggaccttaacctattgtcttgaaatttaggttatctcttagatatcacat	400
gtgattaccccagtgaaacgcgtataagcttacagaaaggaaaaccggttg	450
gctcagtcaaaactggtgcagatttgggctcccctgaatatttgagacat	500
ccctaaaatgaagagatatatacagctaattttgaatgaaaattttaa	550
tcgcaatgaacagtagtagatgagcttttgaagtcctttcaaattatt	600
tggtcttccagttgatattttttattttatataaccagtaccaaataat	650
cttgccatacatttacctttttgagggtgttcaacggaaatccagtgat	700
ttacacattcttggaaacccatcgcttataatacgaactaattttattat	750
gaacaaaggcttggaaaagtatccctactttttacgacgctaaatcatg	800
atacgaaactttaggaagattaacagtcactccataaaatcagaaagtat	850
tcgctaatagtggaaagaaatggttatataa	881

FIG. 51